

ARBOR ASSAYS™
Interactive Assay Solutions™



DetectX®
URINARY RBP
Enzyme Immunoassay Kit

1 Plate Kit Catalog Number KU04-H1

Species Independent

Sample Types Validated:

Human, Rat, Dog and Rhesus Monkey Urine

Please read this insert completely prior to using the product.
For research use only. Not for use in diagnostic procedures.

www.ArborAssays.com   

TABLE OF CONTENTS

Background	3
Assay Principle	4
Related Products	4
Supplied Components	5
Storage Instructions	5
Other Materials Required	6
Precautions	6
Sample Types	7
Sample Preparation	7
Reagent Preparation	8
Assay Protocol	9
Calculation of Results	10
Typical Data	10-11
Validation Data Sensitivity, Linearity, etc.	11-13
Sample Values	14
Warranty & Contact Information	15
Plate Layout Sheet	16

BACKGROUND

Retinol binding protein (RBP) is from a family of structurally related proteins that bind small hydrophobic molecules such as bile pigments, steroids, odorants, etc¹. RBP is a 21 kDa highly conserved, single-chain glycoprotein, consisting of 182 amino acids with 3 disulfide bonds, that has a hydrophobic pocket which binds retinol (vitamin A).

RBP binds retinol in a 1:1 stoichiometry, which serves to not only solubilize retinol but also protect it from oxidation. When in serum, the majority of RBP bound with retinol is reversibly complexed with transthyretin (prealbumin)^{2,3}. This complex then transports retinol to specific receptors of various tissues in the body. Vitamin A status is reflected by serum concentration as it is hemostatically controlled and does not fall until stores are dramatically reduced^{4,5}.

RBP has also been shown to be a useful marker for renal function⁶ as it is totally filtered by the glomeruli and reabsorbed by proximal tubules⁷. This has made urinary RBP (uRBP) a tool to study renal function in heart⁸ or kidney⁹ transplant recipients, type 1 and 2 diabetics¹⁰, and in people exposed to uranium from mining operations¹¹. Measurement of uRBP levels has also been useful in detection and characterization of diseases including hypertension¹² and certain cancers^{13,14}.

1. Blaner WS., "Retinol binding protein: the serum transport protein for vitamin A.", *Endocr Rev.* 1989, 10(3):308-16.
2. Wolf G., "Multiple functions of vitamin A.", *Physiol. Rev.* 1984, 64(3):873-937.
3. Petersen PA., "Characteristics of a Vitamin A-transporting Protein Complex Occuring in Human Serum.", *J. Biol. Chem.* 1971 246:34-43.
4. Goodman DS., Blaner WS., in *The Retinoids*, eds. Sporn MB., Roberts AB., Goodman DS., (Orlando: Academic Press, 1984) vol.2, 1-39.
5. Olson JA., "Vitamin A, retinoids and carotenoids.", *Modern Nutrition in Health and Disease*, eds. Shils ME, Olson JA, Shike M, (Philadelphia: Lea & Febiger, 1994) 8th ed. 287-307.
6. Peterson PA., and Berggård I., "Isolation and Properties of a Human Retinol-transporting Protein.", *J. Biol. Chem.* 1971, 246:25-33.
7. Bernard AM., and Lauwerys RR., "Retinol Bining Protein in Urine: A More Practical Index than Urinary β 2-microglobulin for the Screening of Renal Tubular Function.", *Clin. Chem.* 1981, 27:1781-2.
8. Câmara NO., Matos AC., Rodrigues DA., et al., "Early Detection of Heart Transplant patients with Increased Risk of Cyclosporin Nephrotoxicity.", *The Lancet*, 2001, 357: 856-857.
9. Hosaka B., Park SI., Felipe CR., Garcia RG., Machado PG., Pereira AB., Tedesco-Silva H., and Medina Pestana JO., "Predictive Value of Urinary Retinol Binding Protein for Graft Dysfunction after Kidney Transplantation." *Transplant Proc.*, 2003, 35:1341-1343.
10. Abahusain MA., Wright J., Dickerson JWT., and de Vol EB., "Retinol, α -Tocopherol and Carotenoids in Diabetes." *Am J Clin Nutr.*, 2004, 79:218-25.
11. Wyatt SA., Reitz LV., Croley TR., Hawkins D., Barrett E., McKeown A., Powell N., West A., Hamner T. and Royster MO., "Biological Monitoring of Uranium Exposure in South Central Virginia." *J Exp Sci Environ Epidemi.*, 2008, 18:59-75.
12. Bang LE., Holm J., and Svendsen TL. "Retinol-Binding protein and Transferrin in Urine: New Markers of Renal Function in Essential Hypertension and White Coat Hypertension" *Am. J. Hypertens*, 1996, 9:1024-1028.
13. Moise AR, Noy N, Palczewski K, Blaner WS. "Delivery of retinoid-based therapies to target tissues." *Biochemistry*, 2007, 9:1024-1028.
14. Gavrilov V., Yermiahu T., Gorodischer R., "Renal pathology and retinol status in multiple myeloma patients." *Kidney Int.*, 2006, 69(1):173-7.

ASSAY PRINCIPLE

The DetectX® Urinary Retinol Binding Protein (RBP) kit is designed to quantitatively measure RBP present in urine samples. Please read the complete kit insert before performing this assay. A RBP standard is provided to generate a standard curve for the assay and all samples should be read off the standard curve. Standards or diluted samples are pipetted into a clear microtiter plate coated with an antibody to capture rabbit antibodies. A RBP-peroxidase conjugate is added to the standards and samples in the wells. The binding reaction is initiated by the addition of the RBP polyclonal antibody to each well. After an hour incubation the plate is washed and substrate is added. The substrate reacts with the bound RBP-peroxidase conjugate. After a short incubation, the reaction is stopped and the intensity of the generated color is detected in a microtiter plate reader capable of measuring 450 nm wavelength. The concentration of the RBP in the sample is calculated, after making a suitable correction for the dilution of the sample, using software available with most plate readers.

RELATED PRODUCTS

Kits	Catalog No.
BCA Protein Dual Range Colorimetric Detection Kit	KO41-H1
Human Cystatin C EIA Kit	KO12-H1
Serum Creatinine Detection Kit	KBO2-H1
Thiol Detection Kit	KO05-F1
Urea Nitrogen (BUN) Detection Kits	KO24-H1/H5
Urinary Creatinine Detection Kits	KO02-H1/H5



SUPPLIED COMPONENTS

Coated Clear 96 Well Plate

A clear plastic microplate with break-apart strips coated with goat anti-rabbit IgG.
One Plate Catalog Number X016-1EA

RBP Standard

A stock solution of native human RBP at 20 µg/mL.
60 µL Catalog Number C014-60UL

DetectX® RBP Antibody

A polyclonal antibody specific for RBP.
3 mL Catalog Number C011-3ML

DetectX® RBP-Peroxidase Conjugate

A RBP-peroxidase conjugate.
3 mL Catalog Number C015-3ML

Assay Buffer

28 mL Catalog Number X005-28ML

Wash Buffer Concentrate

A 20X concentrate that should be diluted with deionized or distilled water.
30 mL Catalog Number X007-30ML

TMB Substrate

11 mL Catalog Number X019-11ML

Stop Solution

A 1N hydrochloric acid solution. **CAUSTIC.**
5 mL Catalog Number X020-5ML

Plate Sealer

1 each Catalog Number X002-1EA

STORAGE INSTRUCTIONS

All components of this kit should be stored at 4°C until the expiration date of the kit.

OTHER MATERIALS REQUIRED

Distilled or deionized water.

A microplate shaker and a microplate washer.

Repeater pipet with disposable tips capable of dispensing 25 μ L, 50 μ L, and 100 μ L.

Colorimetric 96 well microplate reader capable of reading optical density at 450 nm.

Software for converting raw relative optical density readings from the plate reader and carrying out four parameter logistic curve (4PLC) fitting. Contact your plate reader manufacturer for details.

PRECAUTIONS

As with all such products, this kit should only be used by qualified personnel who have had laboratory safety instruction. The complete insert should be read and understood before attempting to use the product.

The antibody coated plate needs to be stored desiccated. The silica gel pack included in the foil ziploc bag will keep the plate dry. The silica gel pack will turn from blue to pink if the ziploc has not been closed properly.

The RBP Standard is purified from a human source and as such, should be treated as potentially hazardous. Proper safety procedures must be followed.

This kit utilizes a peroxidase-based readout system. Buffers, including other manufacturers Wash Buffers, containing sodium azide will inhibit color production from the enzyme. Make sure **all** buffers used for samples are **azide free**. Ensure that any plate washing system is rinsed well with deionized water prior to using the supplied Wash Buffer as prepared on Page 8.

The Stop Solution is acid. The solution should not come in contact with skin or eyes. Take appropriate precautions when handling this reagent.



SAMPLE TYPES

This assay has been fully validated for human urine samples and tested in rat, dog and rhesus monkey urines. Samples containing visible particulate should be centrifuged prior to using.

RBP is a highly conserved protein and we have shown that this kit may measure RBP from sources other than human. Please see page 14 for details of other urine samples tested. The end user should evaluate recoveries of RBP in other urine samples being tested.

For measuring retinol binding protein in serum or plasma samples, please refer to the DetectX® Serum RBP Immunoassay kit, Catalog Number K004-H1.

SAMPLE PREPARATION

Samples must be diluted 1:2 by adding one part of urine to one part Assay Buffer prior to running in the kit. Any samples with RBP concentrations greater than the standard curve range should be diluted further with Assay Buffer to obtain readings within the standard curve. Samples that are too dilute to be measured should be concentrated prior to measuring in the assay.

Use all samples within 2 hours of dilution.

REAGENT PREPARATION


Allow the kit reagents to come to room temperature for 30 minutes. Ensure that all samples have reached room temperature and have been diluted as appropriate prior to running them in the kit.

Wash Buffer Preparation

Dilute Wash Buffer Concentrate 1:20 by adding one part of the concentrate to nineteen parts of deionized water. Once diluted this is stable at room temperature for 3 months.

Standard Preparation

Label tubes #1 through #5. Briefly spin vial of standard in a microcentrifuge to ensure contents are at bottom of vial. Pipet 475 μL of Assay Buffer into tube #1 and 300 μL into tubes #2 to #5. Carefully add 25 μL of the RBP stock solution to tube #1 and vortex completely. Take 100 μL of the RBP solution in tube #1 and add it to tube #2 and vortex completely. Repeat these serial dilutions for tubes #3 through #5. The concentration of RBP in tubes 1 through 5 will be 1,000, 250, 62.5, 15.625, and 3.906 ng/mL.



Use all Standards within 2 hours of preparation.

	Std 1	Std 2	Std 3	Std 4	Std 5
Assay Buffer Volume (μL)	475	300	300	300	300
Addition	Stock	Std 1	Std 2	Std 3	Std 4
Volume of Addition (μL)	25	100	100	100	100
Final Conc (ng/mL)	1,000	250	62.5	15.625	3.906

ASSAY PROTOCOL

We recommend that all standards and samples be run in duplicate to allow the end user to accurately determine RBP concentrations.

1. Use the plate layout sheet on the back page of the kit insert to aid in proper sample and standard identification. Determine the number of wells to be used and return unused wells to foil pouch with desiccant. Seal the ziploc plate bag and store at 4°C.
2. Pipet 50 µL of samples or standards into wells in the plate. Pipet 75 µL of Assay Buffer into the non-specific binding (NSB) wells. Pipet 50 µL of Assay Buffer into the maximum binding (BO or Zero standard) wells.
3. Add 25 µL of the DetectX® RBP-peroxidase conjugate to each well, using a repeater pipet.
4. Add 25 µL of the DetectX® RBP Antibody solution to each well, **except the NSB wells**, using a repeater pipet.
5. Gently tap the sides of the plate to ensure adequate mixing of the reagents. Cover the plate with the plate sealer and shake at room temperature for 1 hour.
6. Aspirate the plate and wash each well 4 times with 300 µL wash buffer. Tap the plate dry on clean absorbent towels.
7. Add 100 µL of the TMB Substrate to each well, using a repeater pipet.
8. Incubate the plate at room temperature for 30 minutes without shaking.
9. Add **50 µL** of the Stop Solution to each well, using a repeater pipet.
10. Read the optical density generated from each well in a plate reader capable of reading at 450 nm.
11. Use the plate reader's built-in 4PLC software capabilities to calculate RBP concentration for each sample.

NOTE: If you are using only part of a strip well plate, at the end of the assay throw away the used wells and retain the plate frame for use with the remaining unused wells.

CALCULATION OF RESULTS

Average the duplicate OD readings for each standard and sample. Create a standard curve by reducing the data using the 4PLC fitting routine on the plate reader, after subtracting the mean OD's for the NSB. The sample concentrations obtained, calculated from the %B/B0 curve, should be multiplied by the dilution factor to obtain neat sample values. Sample RBP values should be normalized to creatinine levels by running the same samples in the DetectX® Urinary Creatinine Detection Kit, K002-H1.

Or use the online tool from MyAssays to calculate the data:

www.myassays.com/arbor-assays-retinol-binding-protein-urinary-eia-kit.assay



TYPICAL DATA

Sample	Mean OD	Net OD	% B/B0	RBP Conc. (ng/mL)
NSB	0.098	0	-	-
Standard 1	0.194	0.096	10.7	1,000
Standard 2	0.319	0.221	24.5	250
Standard 3	0.516	0.418	46.4	62.5
Standard 4	0.740	0.642	71.3	15.625
Standard 5	0.907	0.809	89.8	3.906
B0	0.999	0.901	100	0
Sample 1	0.318	0.220	24.4	251.1
Sample 2	0.598	0.500	55.4	38.25

Always run your own standard curve for calculation of results. Do not use this data.

Conversion Factor: 1 ng/mL of human RBP is equivalent to 47.62 pM RBP.



*The MyAssays logo is a registered trademark of MyAssays Ltd.


ARBOR
ASSAYS



WEB INSERT 160902

EXPECT ASSAY ARTISTRY

Typical Standard Curve



Always run your own standard curves for calculation of results. Do not use this data.

VALIDATION DATA

Sensitivity and Limit of Detection

Sensitivity was calculated by comparing the OD's for twenty wells run for each of the B0 and standard #5. The detection limit was determined at two (2) standard deviations from the B0 along the standard curve.

Sensitivity was determined as 2.90 ng/mL.

The Limit of Detection for the assay was determined in a similar manner by comparing the OD's for twenty replicates for each of the zero standard and a low concentration human urine sample.


Limit of Detection was determined as 4.09 ng/mL*.

* Note: Due to the dilute nature of this sample it was run neat instead of being diluted 1:2.

Linearity

Linearity was determined by taking two human urine samples diluted 1:2, one with a low diluted RBP level of 14.5 ng/mL and one with a higher diluted level of 299.2 ng/mL, and mixing them in the ratios given below. The measured concentrations were compared to the expected values based on the ratios used.

High Urine	Low Urine	Expected Conc. (ng/mL)	Observed Conc. (ng/mL)	% Recovery
80%	20%	242.2	227.9	94.1
60%	40%	185.3	180.4	97.3
40%	60%	128.4	120.5	93.9
20%	80%	71.5	70.0	98.0
Mean Recovery				95.8%



Intra Assay Precision

Four human urine samples were diluted 1:2 with Assay Buffer and run in replicates of 8 in an assay. The mean and precision of the calculated RBP concentrations were:

Sample	RBP Conc. (ng/mL)	%CV
1	14.8	4.5
2	28.0	7.5
3	53.5	7.3
4	323.7	2.1

Inter Assay Precision

Four human urine samples were diluted 1:2 with Assay Buffer and run in duplicates in twenty-one assays run over multiple days by three operators. The mean and precision of the calculated RBP concentrations were:

Sample	RBP Conc. (ng/mL)	%CV
1	14.4	9.1
2	28.0	9.2
3	50.8	8.1
4	309.2	14



SAMPLE VALUES

Fourteen random human urine samples were tested in the assay. Values ranged from 6.8 to 788.5 ng/mL with a mean of 114.7 ng/mL. These samples were also run in the DetectX[®] Urinary Creatinine Detection Kit, K002-H1, and the RBP levels normalized to creatinine levels. Normalized values ranged from 35.9 to 573.9 µg RBP/g creatinine.

Normal ranges for urinary RBP are < 130 µg RBP/g creatinine for individuals under 50 years of age and < 172 µg RBP/g creatinine for those equal to or older than 50¹⁵.

Other Species

We have tested a range of urines from other species for RBP. These include rat, dog and rhesus monkey urine. Because of the difficulty in obtaining urine from animals with known medical history, the values we obtained may not be representative of normal or diseased states. Rat urine diluted 1:2 with Assay Buffer had a neat sample value of 43.37 ng/mL and when normalized to urinary creatine gave a reading of 176.1 µg/g creatinine.

Samples of urine from healthy dog and rhesus monkey read below the 3.906 ng/mL standard. We therefore concentrated these samples by freeze drying them and reconstituting them in one-tenth their original volume with Assay Buffer. At that concentration, neat dog urine read at 29.95 ng/mL and when normalized to urinary creatine gave a reading of 32.27 µg/g creatinine. The neat monkey urine read at 10.50 ng/mL and when normalized to urinary creatine gave a reading of 395.5 µg/g creatinine.

15. 84447 Overview: Retinol-Binding Protein, Random, Urine.” Mayo Medical Laboratories: Reference Laboratory services for hospitals worldwide. Web. 02 Sept. 2009.
www.mayomedicallaboratories.com/test-catalog/Overview/84447

LIMITED WARRANTY

Arbor Assays warrants that at the time of shipment this product is free from defects in materials and workmanship. This warranty is in lieu of any other warranty expressed or implied, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose.

We must be notified of any breach of this warranty within 48 hours of receipt of the product. No claim shall be honored if we are not notified within this time period, or if the product has been stored in any way other than outlined in this publication. The sole and exclusive remedy of the customer for any liability based upon this warranty is limited to the replacement of the product, or refund of the invoice price of the goods.

CONTACT INFORMATION

For details concerning this kit or to order any of our products please contact us:

Arbor Assays

1514 Eisenhower Place
Ann Arbor, Michigan 48108 USA

Phone: 734-677-1774

Fax: 734-677-6860

Web: www.ArborAssays.com

E Mail Addresses:

Info@ArborAssays.com

Orders@ArborAssays.com

Technical@ArborAssays.com

Contracts@ArborAssays.com



OFFICIAL SUPPLIER TO ISWE

Arbor Assays and the International Society of Wildlife Endocrinology (ISWE) signed an exclusive agreement for Arbor Assays to supply ISWE members with EIA kits for wildlife conservation research.


DetectX[®], ThioStar[®] and the Arbor Assays logo are all registered trademarks.

www.ArborAssays.com

15

WEB INSERT 160902





 Printed on Forest Stewardship Council certified paper