

Alpha Galactosidase (α -Gal) Activity Assay Kit (Fluorometric)

04/18

(Catalog # K407-100; 100 assays; Store at -20°C)

I. Introduction:

Alpha-Galactosidase (α -Gal; EC 3.2.1.22) hydrolyzes alpha-galactosyl moieties found in glycolipids and glycoproteins. In mammals, α -Gal hydrolyzes poly- and oligosaccharides commonly found in dietary sources that are difficult to digest. Therefore, α -Gal is used in dietary supplements that help to reduce the production of intestinal gases due to consumption of certain foods. It is known total α -Gal activity is due to two major isozymes with unique, yet different thermostability profiles. Alpha-Galactosidase A, is thermolabile and represents approximately 90% of total α -Gal activity. Fabry Disease, a lysosomal disease disorder, is characterized by mutations in alpha-Galactosidase A. These mutations cause abnormal accumulation of glycosphingolipids in lysosomes. BioVision's Alpha Galactosidase Activity Assay Kit provides a simple, rapid way to monitor total α -Gal activity in wide variety of biological samples. In this kit, α -Gal cleaves a synthetic specific substrate releasing a fluorophore, which can be easily quantified (Ex/Em= 360/445 nm). The assay is specific, sensitive and can detect as low as 0.1 μ U of α -Galactosidase activity.



II. Applications:

- Measurement of α -Galactosidase activity in various samples

III. Sample Type:

- Tissue Homogenates: kidney, etc.
- Cell Lysates: U937, etc.
- Biological fluids: Saliva, etc.

IV. Kit Contents:

Components	K407-100	Cap Code	Part Number
α -Gal Assay Buffer	25 ml	NM	K407-100-1
α -Gal Stop Buffer	25 ml	WM	K407-100-2
α -Gal Substrate	220 μ l	Blue	K407-100-3
4-Methylumbelliferone Standard	35 μ l	Yellow	K407-100-4
α -Gal Positive Control	1 vial	Green	K407-100-5

V. User Supplied Reagents and Equipment:

- Multi-well spectrophotometer (ELISA reader)
- 96-well white plate with flat bottom is preferred for this assay. 96-well clear plate can also be used.
- Dounce Tissue Homogenizer (Cat. #1998)

VI. Storage Conditions and Reagent Preparation:

Store kit at -20°C, protect from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay. Upon opening, use within two months.

- **α -Gal Assay Buffer and Stop Buffer:** Store at 4 °C or -20 °C. Bring to 37 °C before use.
- **α -Gal Substrate:** Light sensitive. Thaw at room temperature. Store at -20 °C.
- **4-Methylumbelliferone Standard (5 mM):** Light sensitive. Thaw at room temperature. Store at -20 °C.
- **α -Gal Positive Control:** Reconstitute with 20 μ l α -Gal Assay Buffer and mix thoroughly. Store at -20 °C. Keep on ice while in use. Use within two months.

VII. α -Gal Activity Assay Protocol:

1. **Sample Preparation: For tissue and cells:** Homogenize tissue (10 mg) or pelleted cells ($\sim 5 \times 10^5$) with 100 μ l ice-cold α -Gal Assay Buffer and keep on ice for 10 min. Centrifuge samples at 12,000 x g at 4 °C for 10 min. and collect the supernatant. Dilute the supernatant 10-20 fold in α -Gal Assay Buffer. Add 2-10 μ l of diluted samples into a 96-well plate that will be designated as Sample(s). **For biological fluids:** Undiluted fluids can be added directly to the well. Add 2-10 μ l of samples into well(s) in a 96-well plate that will be designated as Samples. **For Reagent Background Control:** add same volume of α -Gal Assay Buffer in parallel well(s). **For Positive Control:** dilute reconstituted α -Gal Positive Control 1:10 fold with α -Gal Assay Buffer prior to the assay and add 2-6 μ l of diluted α -Gal Positive Control into desired wells(s). Adjust the volume of Positive Control, Sample(s), and Reagent Background Control to **40 μ l/well** with α -Gal Assay Buffer.

Note:

- a. We suggest using 3-5 different volumes of the samples per well to ensure the readings are within the standard curve range and the progress curve rates are within the linear range.
- b. Do not store unused diluted α -Gal Positive Control.

2. **Standard Curve Preparation:** Prepare a 100 μ M 4-Methylumbelliferone (4-MU) Standard by adding 10 μ l of 5 mM 4-MU to 490 μ l α -Gal Assay Buffer in amber tube. Further dilute the 100 μ M Standard solution 5-fold by adding 20 μ l of 100 μ M 4-MU to 80 μ l α -Gal Assay Buffer to generate 20 μ M 4-MU Standard. Add 0, 2, 4, 6, 8, 10 μ l of 20 μ M 4-MU standard into a series of wells to generate 0, 40, 80, 120, 160, 200 pmol/well of 4-MU Standard respectively. Adjust the volume to **60 μ l/well** with α -Gal Assay Buffer.

Note: Equilibrate the α -Gal Assay Buffer to 37 °C prior to the assay.

3. Substrate Hydrolysis: Prepare sufficient volume of 10-fold dilution of α -Gal Substrate (i.e. Dilute 4 μ l of α -Gal stock Substrate with 36 μ l of α -Gal Assay Buffer), vortex briefly. Add 20 μ l of diluted α -Gal Substrate to each well containing the test Sample(s), Positive Control and Reagent Background Control. *The total volume in each well (i.e. Samples, Positive Control and Reagent Background Control) should be 60 μ l.* **Mix well and incubate at 37 °C for 2 hours, avoid light.** After incubation, add 200 μ l of α -Gal Stop Buffer to each well containing Sample(s), Positive Control, Reagent Background Control and Standards. Mix well.

Note:

a. Equilibrate the α -Gal Stop Buffer to 37 °C prior to the assay.

b. Standards can be prepared at the end of the incubation time, and measured in end-point mode.

4. Measurement: Measure fluorescence intensity (Ex/Em= 360/445 nm) at 37°C using an end-point setting.

5. Calculation: Subtract 0 Standard reading from all Standard readings. Plot the 4-MU Standard Curve; subtract the Reagent Background Control reading from all Sample readings. Apply sample Δ RFU to 4-MU Standard Curve to obtain the corresponding pmol of product formed (**B**, in pmol) and calculate the activity of α -Galactosidase activity in the sample as:

$$\text{Specific Sample } \alpha\text{-Galactosidase Activity} = \frac{B}{(2 \times V \times P) \times D} \text{ (pmol/h/mg} \equiv 0.0167 \text{ } \mu\text{U/mg)}$$

Where: **B** = 4-MU amount in sample well from Standard Curve (pmol)

2 = Reaction time (hour)

V = Sample volume added into the reaction well (ml)

P = Initial Sample Concentration in mg-protein/ml (mgP/ml)

D = Sample Dilution Factor

1 pmol/h = 0.0167 pmol/min \equiv 0.0167 μ U

Unit Definition: One unit of α -Galactosidase activity is the amount of enzyme that generates 1.0 μ mol of 4-Methylumbelliferone per min., at pH 4.5 at 37 °C.




Figure: (a) 4-Methylumbelliferone Standard Curve. Results are from multiple experiments. (b) α -Galactosidase Activity in Mouse Kidney Tissue Extracts (1 μ g protein) and U937 Cell Lysates (0.2 μ g protein). (c) Measurement of α -Galactosidase Activity in undiluted Human Pooled Saliva (5 μ l). All assays were performed following kit protocol.

VIII. RELATED PRODUCTS:

β -Galactosidase Activity Assay Kit (K821)
 β -Galactosidase Staining Kit (K802)
 β -Galactosidase Inhibitor Screening Kit (K827)
 α -L-Fucosidase Activity Assay Kit (Fluorometric) (K542)
 α -L-Fucosidase (FUCA1) Assay Kit (Colorimetric)
 β -glucuronidase Activity Assay Kit (K514)
 β -N-Acetylglucosaminidase Activity Assay Kit (Colorimetric) (K733)
 EZClick™ O-GlcNAc Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence) (K714)
 EZClick™ Sialic Acid (ManAz) Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence) (K441)
 EZClick™ O-GalNAc Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence)
 Dounce Tissue Homogenizer (1998)

FOR RESEARCH USE ONLY! Not to be used on humans