

## **HDAC Cell-Based Activity Assay Kit**

Item No. 600150



**Customer Service** 800.364.9897 \* **Technical Support** 888.526.5351

[www.caymanchem.com](http://www.caymanchem.com)

# TABLE OF CONTENTS

GENERAL INFORMATION	3	Materials Supplied
	4	Precautions
	4	If You Have Problems
	4	Storage and Stability
	5	Materials Needed but Not Supplied
INTRODUCTION	6	Background
	6	About This Assay
PRE-ASSAY PREPARATION	7	Cell Culture Preparation
	7	Reagent Preparation
ASSAY PROTOCOL	9	Plate Set Up
	9	Performing the Assay
ANALYSIS	11	Calculations
	13	Performance Characteristics
RESOURCES	15	Troubleshooting
	16	References
	17	Related Products
	18	Warranty and Limitation of Remedy
	19	Plate Template
	20	Notes

# GENERAL INFORMATION

## Materials Supplied

Kit will arrive packaged as a -80°C kit. For best results, remove components and store as stated below.

Item Number	Item	Quantity/Size	Storage
10006389	HDAC Assay Buffer (10X)	1 vial/5 ml	-20°C
600009	Nonidet P-40 Assay Reagent (10%)	1 vial/500 µl	4°C
10006393	HDAC Deacetylated Standard	1 vial/400 µl	-20°C
10011617	HDAC1 Positive Control	1 vial/50 µl	-80°C
600155	Cell-Based Assay HDAC Substrate	1 vial/100 µl	-20°C
10006394	HDAC Developer	1 vial/20 mg	-20°C
10006391	HDAC Trichostatin A	1 vial/250 µl	-20°C
400092	96-Well Clear Bottom Black Culture Plate	1 plate	Room Temperature

If any of the items listed above are damaged or missing, please contact our Customer Service department at (800) 364-9897 or (734) 971-3335. We cannot accept any returns without prior authorization.

*Sold under license from CycLex® Co., Ltd (Patent Nos. EP 1243658, US 7033778, US 7256013, and JP 4267043).*



**WARNING:** This product is for laboratory research use only: not for administration to humans. Not for human or veterinary diagnostic or therapeutic use.

## Precautions

**Please read these instructions carefully before beginning this assay.**

**For research use only. Not for human or diagnostic use.**

## If You Have Problems

### Technical Service Contact Information

**Phone:** 888-526-5351 (USA and Canada only) or 734-975-3888

**Fax:** 734-971-3641

**Email:** techserv@caymanchem.com

**Hours:** M-F 8:00 AM to 5:30 PM EST

In order for our staff to assist you quickly and efficiently, please be ready to supply the lot number of the kit (found on the outside of the box).

## Storage and Stability

The kit will perform as specified if components are stored as directed on page 3 and used before the expiration date indicated on the outside of the box.

## Materials Needed But Not Supplied

1. Cell line such as HeLa cells or MCF-7 cells (can be obtained from ATCC); other cell lines may also be used
2. A fluorometer with the capacity to measure fluorescence using an excitation wavelength of 340-360 nm and an emission wavelength of 440-465 nm
3. Adjustable pipettes and a repeating pipettor
4. A source of pure water; glass distilled water or HPLC-grade water is acceptable

## Background

Nucleosomes, which fold chromosomal DNA, contain two molecules each of the core histones H2A, H2B, H3, and H4. Almost two turns of DNA are wrapped around this octameric core, which represses transcription.<sup>1</sup> The histone amino termini extend from the core, where they can be modified post-translationally by acetylation, phosphorylation, ubiquitination, and methylation, affecting their charge and function.

Acetylation of the  $\epsilon$ -amino groups of specific histone lysines is catalyzed by histone acetyltransferases (HATs) and correlates with an open chromatin structure and gene activation. Histone deacetylases (HDACs) catalyze the hydrolytic removal of acetyl groups from histone lysine residues and correlates with chromatin condensation and transcriptional repression.<sup>2,3</sup> Therefore, HDAC inhibition results in transcriptional activation through the conformational relaxation of DNA. Changes in the transcription of key genes has linked HDAC inhibitors to blocking angiogenesis and cell cycling, and promoting apoptosis and differentiation. By targeting these key components of tumor proliferation, HDAC inhibitors are currently being explored as potential anticancer agents.<sup>4,6</sup>

## About This Assay

Cayman's HDAC Cell-Based Assay Kit provides an easy tool for studying HDAC activity modulators in whole cells. By using a cell-permeable HDAC substrate, the activity of various protein lysine-specific deacetylases including HDAC1-containing complexes can be measured in intact cells in a simple and homogenous manner. The fluorescence of the deacetylated reaction product can be analyzed using a plate reader or a fluorometer with excitation wavelengths of 340-360 nm and emission wavelengths of 440-465 nm. An HDAC inhibitor, trichostatin A, is included for checking specificity of the HDAC reaction. This assay parallels Cayman's HDAC Activity Assay Kit (Item No. 10011563), which uses a nuclear extract rather than whole cells for the assay. Together, both assays will help to identify whether an inhibitor/activator has a direct effect on the enzyme.

## Cell Culture Preparation

1. Seed cells in the provided clear bottom black 96-well plate at a density of  $(2 \times 10^4)$  -  $(5 \times 10^4)$  cells/well in 100  $\mu$ l of culture medium. For a blank control, add 100  $\mu$ l of culture medium without cells to 2-3 wells of the plate. *NOTE: If you run a standard curve with the assay, add 100  $\mu$ l of culture medium without cells to 12 wells of the plate.* Culture the cells overnight or until the cells reach 80% confluence.
2. Treat the cells with or without compounds to be tested. We recommend that each treatment be performed in duplicate.
3. Following addition of test compounds, continue to culture the cells in a CO<sub>2</sub> incubator at 37°C for 24-48 hours, or for a period of time according to your typical experimental protocol.

## Reagent Preparation

### 1. HDAC Assay Buffer (10X)

Dilute 5 ml of HDAC Assay Buffer (10X) (Item No. 10006389) with 45 ml of UltraPure water. This final Assay Buffer should be used for diluting the HDAC Deacetylated Standard (Item No. 10006393), HDAC1 Positive Control (Item No. 10011617), Cell-Based Assay HDAC Substrate (Item No. 600155), and for dissolving the HDAC Developer (Item No. 10006394). The diluted buffer is stable for six months at -20°C.

## 2. HDAC Deacetylated Standard

The HDAC Deacetylated Standard (Item No. 10006393) vial contains 400 µl of 2.1 mM deacetylated standard in dimethylsulfoxide (DMSO). The deacetylated standard is used to prepare a standard curve for quantitative determination of HDAC activity.

To run a standard curve, obtain six clean test tubes and label them #1 through #6. Add 1.4 ml of Assay Buffer into tube #1 and 500 µl into tubes #2-6. Transfer 20 µl of the HDAC Deacetylated Standard to tube #1 and mix thoroughly. The concentration of this Standard, the first point on the standard curve, is 30 µM. Serially dilute the Standard by removing 500 µl from tube #1 and placing it into tube #2; mix thoroughly. Next remove 500 µl from tube #2 and place it into tube #3; mix thoroughly. Repeat the procedure for tubes #4-5. The concentration of these standards will be 30, 15, 7.5, 3.75, and 1.88 µM, respectively. Do not add any standard to tube #6. This tube will be your blank.

## 3. HDAC1 Positive Control

The HDAC1 Positive Control (Item No. 10011617) vial contains 50 µl of human recombinant HDAC1. Dilute 10 µl of the control with 190 µl of diluted Assay Buffer. The diluted HDAC1 is stable for four hours when stored on ice.

## 4. HDAC Trichostatin A

The HDAC Trichostatin A (Item No. 10006391) vial contains 250 µl of 0.21 mM Trichostatin A. Trichostatin A is an HDAC inhibitor. Dilute 50 µl of Trichostatin A stock with 450 µl of diluted Assay Buffer. A 10 µl aliquot in the assay results in a final concentration of 2.1 µM. At this concentration, HDAC activity will be completely inhibited.

## 5. Cell-Based Assay HDAC Substrate

The Cell-Based Assay HDAC Substrate (Item No. 600155) vial contains 100 µl of concentrated Boc-Lys(AC)-AMC in DMSO. Prior to assaying, dilute 100 µl of the Substrate with 1 ml of diluted Assay Buffer.

## 6. HDAC Lysis/Developer Mixture

*NOTE: this HDAC Lysis/Developer Mixture should be prepared immediately before use.*

To one vial of Developer, add 4.5 ml of diluted Assay Buffer, 72 µl of HDAC Trichostatin A, and 500 µl of Nonidet P-40 Assay Reagent (10%) (Item No. 600009).

## ASSAY PROTOCOL

### Plate Set Up

There is no specific pattern for using the wells on the plate. However, a positive control in duplicate has to be assayed with the sample. We suggest that each sample be assayed at least in duplicate and to have two wells designated as background wells. We also recommend assaying each sample in the presence and absence of the HDAC inhibitor to allow for the correction of HDAC-independent fluorescence. Record the contents of each well on the template sheet provided on page 19.

### Performing the Assay

#### Pipetting Hints

- Use different tips to pipette each reagent.
- Before pipetting each reagent, equilibrate the pipette tip in that reagent (*i.e.*, slowly fill the tip and gently expel the contents, repeat several times).
- Do not expose the pipette tip to the reagent(s) already in the well.

## Procedure

1. Centrifuge the plate in a plate centrifuge at 500 x g for five minutes.
2. Aspirate the culture medium.
3. Add 200 µl of diluted Assay Buffer to each well and centrifuge the plate at 500 x g for five minutes.
4. Aspirate the supernatant.
5. Add 90 µl of culture medium or positive control to non-inhibited sample wells. Add 80 µl of culture medium plus 10 µl of Trichostatin A to appropriate control wells to test assay specificity. Initiate the HDAC reactions by adding 10 µl of diluted HDAC Substrate to each well.
6. Incubate the plate at 37°C for two hours for optimal development (can be 1-3 hours according to your schedule).
7. Add 100 µl standards prepared above to appropriate wells. Skip this step if you are not running a standard curve with the assay.
8. Add 50 µl of the Lysis/Developer Solution to each well.
9. Shake the plate on a plate shaker for a 1-2 minutes.
10. Incubate the plate for 15 minutes at 37°C.
11. Read the fluorescent intensity of each well (excitation = 340-360 nm; emission = 440-460 nm).

## ANALYSIS

### Calculations

#### Plot the Standard Curve

1. Determine the average fluorescence of the standards. Subtract the fluorescence value of the blank (standard tube #6) from itself and all other standards. This is the corrected fluorescence.
2. Plot the corrected fluorescence values (from step 1 above) of each standard as a function of the final concentration of Deacetylated Standards. See Figure 1, on page 12, for a typical standard curve.


#### Determination of HDAC activity

1. Determine the average fluorescence of each sample and sample plus Trichostatin A.
2. Subtract the Trichostatin A sample fluorescence from the non-Trichostatin A sample fluorescence to yield the corrected sample fluorescence.
3. Calculate the HDAC activity using the equation obtained from the linear regression of the standard curve, substituting corrected fluorescence values for each sample.

HDAC Activity (nmol/min/ml) =


$$\frac{[\text{Corrected Sample Fluorescence} - (\text{y-intercept})/\text{slope}] \times \text{dilution}}{15 \text{ minutes}}$$

If you anticipate a high activity of HDAC in the samples, dilution may be required to obtain values that fall on the standard curve.




**Figure 1. Deacetylated substrate standard curve**

## Performance Characteristics



**Figure 2. Inhibition of HDAC activity by SAHA.** MCF-7 cells were seeded in a 96-well clear bottom black plate at a density of  $4 \times 10^4$  cells/well. The next day, cells were treated with different concentrations of SAHA as indicated above. On the third day, cells were processed for the measurement of HDAC activity according to the protocol described in this booklet.



**Figure 3. Effect of substrate incubation time on HDAC activity in HeLa cells and MCF-7 cells.** HeLa cells and MCF-7 cells were seeded in a 96-well clear bottom black plate at a density of  $4 \times 10^4$  cells/well. The next day, cells were processed for the measurement of HDAC activity according to the protocol described in this booklet. One hour of substrate incubation generates measurable fluorescence intensity whereas two hours of substrate incubation time doubles the amount of fluorescence activity indicating linearity of the assay over this time period. Addition of TSA almost completely blocks fluorescence, indicating specificity of the HDAC assay.

## RESOURCES

### Troubleshooting

Problem	Possible Causes	Recommended Solutions
Erratic values; dispersions of duplicates/triplicates	A. Poor pipetting/technique B. Bubble in the well(s)	A. Be careful not to splash the contents of the wells B. Carefully tap the side of the plate with your finger to remove bubbles
Erratic response curve of compound treatments	Unequal number of cells in each well	Make sure each well contains the same number of cells
The fluorometer exhibited 'MAX' values for the wells	The GAIN setting is too high	Reduce the GAIN and re-read
High reading in all wells	Cell density is too high	Plate cells more sparsely



## References

1. Grunstein, M. Histone acetylation in chromatin structure and transcription. *Nature* **389**, 349-352 (1997).
2. Strahl, B.D. and Allis, D. The language of covalent histone modifications. *Nature* **403**, 41-45 (2000).
3. Cheung, W.L., Briggs, D.B., and Allis, C.D. Acetylation and chromosomal functions. *Curr. Opin. Cell Biol.* **12**, 326-333 (2000).
4. Kwon, H.J., Kim, M.S., Kim, M.J., *et al.* Histone deacetylase inhibitor FK228 inhibits tumor angiogenesis. *Int. J. Cancer* **97**, 290-296 (2002).
5. Suenaga, M., Soda, H., Oka, M., *et al.* Histone deacetylase inhibitors suppress telomerase reverse transcriptase mRNA expression in prostate cancer cells. *Int. J. Cancer* **97**, 621-625 (2002).
6. Johnstone, R.W. Histone-deacetylase inhibitors: Novel drugs for the treatment of cancer. *Nat. Rev. Drug Discov.* **1**, 287-299 (2002).

## Related Products

CAY10398 - Item No. 89740  
CAY10433 - Item No. 10005019  
CAY10591 - Item No. 10009797  
HAT Inhibitor Screening Assay Kit - Item No. 10006515  
HDAC Fluorometric Activity Assay Kit - Item No. 10011563  
HDAC1 Inhibitor Screening Assay Kit - Item No. 10011564  
HDAC8 Inhibitor Screening Assay Kit - Item No. 700230  
NCOR2/SMRT (human recombinant) - Item No. 11633  
PAD4 Inhibitor Screening Assay Kit - Item No. 700560  
SIRT1 Direct Fluorescent Screening Assay Kit - Item No. 10010401  
SIRT1 FRET-Based Screening Assay Kit - Item No. 10010991  
SIRT3 Direct Fluorescent Screening Assay Kit - Item No. 10011566  
Splitomicin - Item No. 13168  
Trichostatin A - Item No. 89730

Warranty and Limitation of Remedy

Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman **warrants only** to the original customer that the material will meet our specifications at the time of delivery. Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, indirect, incidental or consequential damages, even if Cayman is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees.

Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a refund of the purchase price, or at Cayman's option, the replacement, at no cost to Buyer, of all material that does not meet our specifications.

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H

## NOTES

This document is copyrighted. All rights are reserved. This document may not, in whole or part, be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form without prior consent, in writing, from Cayman Chemical Company.

©04/02/2015, Cayman Chemical Company, Ann Arbor, MI, All rights reserved. Printed in U.S.A.