

Yeast Cell Lysis Preparation Kit

2015-200

10 preps

2015-400

25 preps

2015-600

100 preps

- For Use with the G NOME[®] DNA Isolation Kit and the RNaid[®] PLUS Kit
- Suitable for Genomic and Plasmid DNA Isolation

Shipping & Storage:

Store all components at room temperature. Once **Spheroplasting Enzyme Mixx** is made store at -70°C .

Revision #: 2015-999-6A01P



*"Check with BIO 101...
We have you time in mind!"*

1070 Joshua Way, Vista CA 92083, USA

Tel: (760) 598-7299 • Fax (760)598-0116

E-mail: technical@bio101.com

Website: www.bio101.com

For technical support call (800) 424-6101

Yeast Cell Lysis Protocol

- Kit Components 3**
 - 2015-200 (10 preps)3
 - 2015-400 (25 preps)3
 - 2015-600 (100 preps)3
- Introduction4**
- Procedure5**
 - How to Monitor Yeast Cell Lysis to Indicate Spheroplast Formation5
- Product Use Limitation & Warranty 8**
 - General Information8

Kit Components

2015-200 (10 preps)

Name	Volume	Catalog #
Yeast Suspension Buffer	10 ml	2015-201
Yeast Enzyme Enhancer	6 ml	2015-202
Spheroplast Indicating Solution	400 μ l	2015-203
Spheroplast Control Solution	400 μ l	2015-204
Yeast Enzyme Salts	120 μ l	2015-205
Spheroplasting Enzyme Solution	800 μ l	2015-206
Spheroplasting Enzyme Mixx	5.7 mg	2015-207

2015-400 (25 preps)

Name	Volume	Catalog #
Yeast Suspension Buffer	25 ml	2015-401
Yeast Enzyme Enhancer	14 ml	2015-402
Spheroplast Indicating Solution	1 ml	2015-403
Spheroplast Control Solution	1 ml	2015-404
Yeast Enzyme Salts	300 μ l	2015-405
Spheroplasting Enzyme Solution	2 ml	2015-406
Spheroplasting Enzyme Mixx	14.25 mg	2015-407

2015-600 (100 preps)

Name	Volume	Catalog #
Yeast Suspension Buffer	100 ml	2015-601
Yeast Enzyme Enhancer	56 ml	2015-602
Spheroplast Indicating Solution	4 ml	2015-603
Spheroplast Control Solution	4 ml	2015-604
Yeast Enzyme Salts	1.2 ml	2015-605
Spheroplasting Enzyme Solution	4 x 2 ml	2015-606
Spheroplasting Enzyme Mixx	4 x 14.25 mg	2015-607

Introduction

The **Yeast Cell Lysis Preparation** Kit contains the necessary reagents for removal of the yeast cell wall in preparation for efficient lysis with the **G NOME** DNA isolation kit (Catalog #'s: 2010-200, 2010-400, 2010-600). The rapid spheroplasting method employed in this protocol is suitable for genomic and plasmid DNA isolation. Transformation and tissue culture with these spheroplasts should not be performed following this method (for this protocol see **Yeast Spheroplast Transformation** Kit catalog #2210-200).

Procedure

1. Start with 100 ml of yeast grown to a density of 2×10^7 cells/ml.
2. Pellet at 600 xg for 5 minutes.
3. Discard supernatant and resuspend pellet in 1 ml of **Yeast Suspension Buffer**. Transfer to a microcentrifuge tube and centrifuge for 10 seconds.
4. Discard supernatant and resuspend in 500 μ l of **Yeast Enzyme Enhancer**.
5. Add 10 μ l **Yeast Enzyme Salts** and 80 μ l **Spheroplasting Enzyme Mixx**.

Before first use combine Spheroplasting Enzyme Solution with Spheroplasting Enzyme Mixx. Mix at room temperature for 5 minutes, aliquot into 80 μ l fractions (the solution often appears 'grainy' at this point). Quick freeze in liquid nitrogen. Store at -70°C until ready to use, only thaw once.

6. Incubate at 37°C until spheroplast formation is complete (usually 12-30 minutes).

How to Monitor Yeast Cell Lysis to Indicate Spheroplast Formation

Start monitoring spheroplast formation after 12 minutes of incubation. Allowing the reaction to proceed after spheroplasting is complete is not beneficial to subsequent procedures.

Place 20 μ l of **Spheroplast Indicating Solution** at one end of a glass slide and 20 μ l of **Spheroplast Control Solution** at the other end. After the cells have incubated for 12 minutes mix 2 μ l of yeast cells with each solution, cover with cover slips, and observe under the microscope. Spheroplasting is sufficient when less than 5% of the control number of cells remain intact in the visual field with the indicating solution relative to the control solution.

At this point cells are easily lysed using the general G NOME protocol.

Notes

Notes

Product Use Limitation & Warranty

Unless otherwise indicated, this product is for research use only. Purchase of BIO 101 products does not grant rights to reproduce, modify, repackage the products or any derivative thereof to third parties. BIO 101 makes no warranty of any kind, expressed or implied, including merchantability or fitness for any particular purpose, except that the products sold will meet our specifications at the time of delivery. Buyer's exclusive remedy and BIO 101's sole liability hereunder shall be limited to, at our option, product credits, refund of the purchase price of, or the replacement of all material(s) that does not meet our specification. By acceptance of the product, Buyer indemnifies and holds BIO 101 harmless against, and assumes all liability for the consequence of its use or misuse by the Buyer, its employees or others, including, but not limited to, the cost of handling. Said refunds or replacement is conditioned on Buyer notifying BIO 101 within thirty (30) days of receipt of product. Failure of Buyer to give said notice within thirty (days) shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

General Information

BIO 101 is a pioneer in developing kits for molecular biology research. We introduced the **GENECLEAN**[®] Kits in 1986 and have since been manufacturing products to bring convenience into your research. Our goal is to make your life easier by simplifying the complexities of lab work.

Technical Support and Ordering Information

(800) 424-6101 (Toll Free - Continental USA)

(760) 598-7299 (Outside USA)

(760) 598-0116 (Fax line)

technical@bio101.com

order@bio101.com

Office Hours:

6:30 am - 6:00 pm P.S.T. (Mon - Fri)

Mailing Address:

BIO 101, Inc.

P.O. Box 2284

La Jolla, CA 92038-2284 USA

Call (800) 424-6101 for Technical Support