

PCR Cycler Check™

For conventional PCR block cyclers

INSTRUCTIONS FOR USE

FOR USE IN RESEARCH AND QUALITY CONTROL

Symbols



Lot No.



Order No.



Expiry date



Storage temperature



Number of reactions



Manufacturer

INDICATION

False negative PCR results or unspecific amplifications might be caused by a PCR cycler defect. Such cases are critical but can be identified by assessing the temperature accuracy of the PCR cycler. However, temperature assessment of a PCR cycler needs special and therefore expensive equipment, such as temperature sensors that measure the temperature homogeneity in a cycler block.

The PCR Cycler Check kit is specifically designed for verifying conventional PCR cyclers, particularly for installation qualification (IQ), operational qualification (OQ) and performance qualification (PQ) as required in various international norms, such as EN ISO 17025, EN 45001, EN ISO 13485, ISO/TS 20836:2007, GLP, GMP, and others.

TEST PRINCIPLE

The PCR Cycler Check kit is based on a temperature sensitive PCR assay to monitor an upper and lower temperature range in one run. The primer sequences in combination with a regular PCR protocol were designed to react extremely sensitive to temperature deviations, temperature homogeneity, precision and timing. Amplification will be altered and indicated with different band pattern at temperature differences of more than 2 °C. The cycler performance is tested at standard PCR settings to represent common applications.

In addition, the pre-adjusted target concentrations are only amplified at high PCR efficiencies as an additional indicator for accurate temperature control of the thermal cycler.

REAGENTS

Each kit contains all reagents required to run the PCR. The expiry date of the unopened package is marked on the package label. The kit components must be stored until use at +2 to +8 °C.

Kit component	Quantity	
	Advance Cat. No. 57-2102	OneStep Cat. No. 57-2103
Validation Reactions	6 strips, 8 vials each, freeze-dried, pre-cast	4 vials, for 25 reactions each, freeze-dried
Caps	6 cap strips, domed	n.a.
Rehydration Buffer	1 vial (1.6 ml)	2 vials (1.6 ml each)
Marker	1 vial (50 µl)	2 vials (50 µl)

The lot specific Certificate of Analysis can be downloaded from our website (www.minerva-biolabs.com).

USER-SUPPLIED CONSUMABLES AND EQUIPMENT

The PCR Cyclor Check kit contains reagents and consumables to perform the cyclor check. Additional consumables and equipment are supplied by the user:

- PCR device for 0.2 ml PCR tubes (relevant only for Cat. No. 57-2102).
- Suitable PCR reaction tubes (relevant only for Cat. No. 57-2103)
- 96-well rack for 0.2 ml PCR tubes (relevant only for Cat. No. 57-2102)
- Microcentrifuge for 8-strips (relevant only for Cat. No. 57-2102) and 2 ml reaction tubes
- Vortexer
- Pipettes with corresponding filter tips

PRECAUTIONS

The PCR Cycler Check™ kit is for in vitro use only. The kit should be used by trained laboratory staff only.

The PCR Cycler Check™ kit does not contain hazardous substances. Remnants can be discarded according to local regulations.

ADDITIONAL NOTES

- ⇒ These instructions must be understood to successfully use the PCR Cycler Check™ kit. The reagents supplied should not be mixed with reagents from different batches but used as an integral unit. The reagents of the kit must not be used beyond their shelf life.
- ⇒ Follow the exact protocol. Any deviation may affect the test method and can affect the results.
- ⇒ Additional control samples are not required. The kit already contains all necessary controls.

PROCEDURE

1A. Reagent preparation for Advance format (Cat. No. 57-2102)


1. Spin down the validation strips to collect the lyophilized material at the bottom of the tube and place the strips in a 96-well rack. Spin down the rehydration buffer.
 2. Carefully remove the protective seal from the validation strips.
 3. Aliquot 25 μ l Rehydration Buffer into each PCR reaction tube. Close the tubes with the provided cap strips.
 4. Incubate for 5 min at room temperature.
 5. Vortex briefly and spin down for 5 sec. Proceed immediately with the PCR.
-

1B. Reagent preparation for OneStep format (Cat. No. 57-2103)

1. Spin down the Validation Tubes and the Rehydration Buffer.
 2. Add 650 μ l of the Rehydration Buffer (blue cap) to the Validation Tube (red cap).
 3. Incubate for 5 min at room temperature.
 4. Vortex briefly and spin down for 5 sec.
 5. Note: Proceed immediately to step 6 or store the mix at < -18 °C. Repeated freezing and thawing must be avoided. We recommend storing the mix in aliquots.
 6. Aliquot 25 μ l of the rehydrated Validation Reagent into each PCR tube.
 7. Close the PCR tubes and spin down briefly. Proceed immediately with the PCR.
-

2. Perform the PCR cyclers test

Place the PCR tubes in the cycler. We recommend the following scheme depending on the cycler block format:



Program the cycler as follows:

Step 1 (pre-incubation): 94 °C for 2 min

Step 2 (amplification):

Cycles 35

Denaturation 94 °C for 30 sec

Annealing T_a for 30 sec (Annealing Temperature (T_a) is provided on the Certificate of Analysis (CoA))

Elongation 72 °C for 30 sec

Step 3:

Hold 4 °C to 8 °C

3. Analysis

1. Prepare a 1.5 % agarose gel, approx. 5 mm thick, with a 5 mm comb.

Load 5 μ l of each PCR reaction. Load 5 μ l of the provided marker in each lane.

2. Note: Loading buffer with dye is already included in the mixes. Thus additional loading buffer or dye is not required.

3. Perform the gel electrophoresis (e.g. 20 min at 100 V).

4. Visualize the PCR results on a UV transilluminator

RESULT INTERPRETATION

The cyclers passed the test if a single band is visible (Fig. 1). The test run is valid but the cycler does not comply with the expected specifications if either no band or two bands are visible.

If no band is visible in any reaction, the experiment should be repeated to exclude a setup mistake. For the re-test, the annealing temperature (T_a) should be reduced by 3 °C to enhance amplification. If the re-test does not show amplification products and the cycler is already suspected to work out of specification, the device should be sent in for service.


If two bands are visible, either the setup of the test was not correct or the cycler is out of specification and should be sent in for service.

Please note, that all PCR reactions must show a uniform result. If this is not the case, most likely one or even more of the Peltier elements have a malfunction. In this case the experiment should be repeated with an adopted loading scheme.

Fragment size	Interpretation
144 bp and 210 bp	annealing temperature too low denaturation temperature ok
144 bp	Cycler test passed successfully
no bands	annealing temperature too high (s. explanation above) or/and denaturation temperature failure

Fig. 1: Gel figure showing results obtained at different annealing temperatures

- A: Marker
- B: Temperature too low
- C: Temperature correct
- D: Temperature too high



APPENDIX

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising from the use, the results of use, or the inability to use this product.

Trademarks

Venor, Onar and AquaScreen are registered trademarks and PCR Cycler Check, ConviFlex, SwabUp, Mycoplasma Off, Food Control, Meat ID, Vegan Control, ExtractNow, WaterShield, Lab Clean and PCR Clean are trademarks of Minerva Biolabs GmbH, Germany.

Related Products

qPCR Cyclor Check™

57-2201 qPCR Cyclor Validation 100 reactions

ConviFlex™ DNAmix

191-025/100/250 PCR Mix with Taq polymerase for conventional and qPCR 25/100/250 reactions

SwabUp™ Lab Monitoring Kits

181-0010/0050 Sample collection and DNA extraction 10/50 samples

182-0010/0050 Sample collection, DNA extraction and PCR system 10/50 samples

Food and Water Assays

11-02-XX-025 Food Control™ qPCR 25 reactions

12-01-005/-020/-040 Meat ID™ Screen 5/20/40 tests

12-02-025/-100 Meat ID™ Halal 25/100 reactions

12-05-025/-100 Vegan Control™ 25/100 reactions

34-2025/-2100/-2250 AquaScreen® qPCR 25/100/250 reactions

Contamination Control Kits for conventional PCR

11-1025/1050/1100/1250 Venor®GeM Classic Mycoplasma Detection Kit 25/50/100/250 tests

11-7024/7048/7096/7240 Venor®GeM Advance Mycoplasma Detection Kit 24/48/96/240 tests

11-8025/8050/8100/8250 Venor®GeM OneStep Mycoplasma Detection Kit 25/50/100/250 tests

12-1025/1050/1100/1250 Onar® Bacteria Detection Kit 25/50/100/250 tests

Contamination Control Kits for qPCR

11-9025/9100/9250 Venor®GeM qEP Mycoplasma Detection Kit 25/100/250 tests

Nucleic Acid Extraction

601-1010/-1050 ExtractNow™ DNA Mini Kit 10/50 extractions

602-1010/-1050 ExtractNow™ Blood DNA Mini Kit 10/50 extractions

603-1010/-1050 ExtractNow™ RNA Mini Kit 10/50 extractions

604-1010/-1050 ExtractNow™ CleanUp Kit 10/50 extractions

605-1010/-1050 ExtractNow™ Plasmid Mini Kit 10/50 extractions

606-1010/-1050 ExtractNow™ Virus DNA/RNA Kit 10/50 extractions

MB Taq DNA Polymerase

53-0050/0100/0200/0250 MB Taq DNA Polymerase (5 U/μl) 50/100/200/250 units

53-1050/1100/1200/1250 MB Taq DNA Polymerase (1 U/μl) 50/100/200/250 units

PCR Clean™

15-2025 DNA Decontamination Reagent, spray bottle 250 ml

15-2200 DNA Decontamination Reagent, refill bottles 4 x 500 ml

15-2201 Wipes 120 wipes in a dispenser box

15-2202 Wipes, refill packs 5 x 120 wipes in a bag

15-2203 Wipes, single wrapped 30 Sachets

Lab Clean™

15-4100 Molecular microbiology lab cleaner, bottled 1 Liter

WaterShield™

15-3015/3020/3050 Water Disinfection Additive for incubators and water baths, 200x concentrate 30 x 5 ml/3 x 50 ml/500 ml

Minerva Biolabs GmbH
Schkopauer Ring 13
D-12681 Berlin, Germany

www.minerva-biolabs.com
Ordering: order@minerva-biolabs.com
Support: support@minerva-biolabs.com

USA & Canada

Minerva Biolabs Inc.
1 Jill Ct., Building 16, Unit 10
Hillsborough, NJ 08844
USA

www.minervabiolabs.us
Ordering: order@minervabiolabs.us
Support: help@minervabiolabs.us

Made in Germany

© 2019 Minerva Biolabs
HB18.05EN