

DATA EXAMPLES

Differentiation of naïve CD4⁺ T cells into Th1 cells is confirmed by intracellular staining for IFN- γ (Figure 1) and secretion of IFN- γ (Figure 2). The corresponding tests for IL-4 (Th2 cell marker) and IL-17 (Th17 cell marker) are negative.




Figure 1: Differentiated Human CD4⁺ T Cells Stain for IFN- γ . Flow cytometry data showing human peripheral blood naïve CD4⁺ T cells without (A, C) and with (B, D) a 5 day differentiation using reagents included in the Human Th1 Cell Differentiation Kit. (A, B) The cells were stained with a Human IFN- γ APC Monoclonal Antibody and a Human IL-4 PE Monoclonal Antibody. (C, D) The cells were stained with a Human IFN- γ APC Monoclonal Antibody and a Human IL-17 PerCP Monoclonal Antibody. Control cultures were used to place the quadrants. All R&D Systems antibodies and corresponding catalog numbers used in this figure are shown below.

SUGGESTED REAGENTS FOR FLOW CYTOMETRY

CATALOG #	DESCRIPTION
IC285A, and IC0041A	Human IFN- γ APC MAb (Clone 25723), Mouse IgG _{2b} , and Mouse IgG _{2b} APC Isotype Control (Clone 133303)
IC204P, and IC002P	Human IL-4 Phycoerythrin MAb (Clone 3007), Mouse IgG ₁ , and Mouse IgG ₁ Phycoerythrin Isotype Control (Clone 11711)
IC3171C, and IC002C	Human IL-17 PerCP MAb (Clone 41802), Mouse IgG ₁ , and Mouse IgG ₁ PerCP Isotype Control (Clone 11711)
FAB3791F, and IC003F	Human CD4 Fluorescein MAb (Clone 11830), Mouse IgG _{2a} , and Mouse IgG _{2a} Fluorescein Isotype Control (Clone 20102)
FC004	Flow Cytometry Fixation Buffer (1X)
FC005	Flow Cytometry Permeabilization/Wash Buffer (1X)

DATA EXAMPLES CONTINUED




Figure 2: Differentiated Human CD4⁺ T Cells Secrete IFN- γ .

Human peripheral blood naïve CD4⁺ T cells were cultured for 5 days using reagents included in the Human Th1 Cell Differentiation Kit. Cytokine expression was determined using the Human IFN- γ Quantikine ELISA Kit, the Human IL-4 Quantikine ELISA kit, and the Human IL-17 Quantikine ELISA Kit. All relevant R&D Systems Quantikine ELISA kits and corresponding catalog numbers are listed below.

SUGGESTED REAGENTS FOR ELISA

CATALOG #	DESCRIPTION
DIF50, or DY285	Human IFN- γ Quantikine ELISA Kit, or Human IFN- γ DuoSet [®]
D4050, or DY204	Human IL-4 Quantikine ELISA Kit, or Human IL-4 DuoSet
D1700, or DY317	Human IL-17 Quantikine ELISA Kit, or Human IL-17 DuoSet

REFERENCES

- Zhu, J. and W.E. Paul (2010) Immunol. Rev. **238**:247.
- Dardalhon, V. *et al.* (2008) J. Autoimmun. **31**:252.
- Seder, R.A. and W.E. Paul (1994) Annu. Rev. Immunol. **12**:635.

CellXVivo™

Human Th1 Cell Differentiation Kit

Catalog Number: CDK001

BACKGROUND

T helper type 1 (Th1) cells are a lineage of CD4⁺ effector T cells that promote cell-mediated immune responses against intracellular viral and bacterial pathogens (1). Enhanced Th1 responses are associated with autoimmune diseases including rheumatoid arthritis, multiple sclerosis, inflammatory bowel disease, and type-1 diabetes (2). Differentiation of CD4⁺ effector T cells into the Th1 lineage is promoted by cytokines such as IL-12 and IFN- γ (3). Th1 cells secrete IFN- γ , IL-10, and TNF- α . The Human Th1 Differentiation Kit contains all necessary components to differentiate human naïve CD4⁺ T cells into Th1 polarized cells.

This package insert must be read in its entirety before using this product. For research use only. Not for use in diagnostic procedures.

MANUFACTURED AND DISTRIBUTED BY:

USA & Canada | R&D Systems, Inc.

614 McKinley Place NE, Minneapolis, MN 55413, USA
TEL: (800) 343-7475 (612) 379-2956 FAX: (612) 656-4400
E-MAIL: info@RnDSystems.com

DISTRIBUTED BY:

UK & Europe | R&D Systems Europe, Ltd.

19 Barton Lane, Abingdon Science Park, Abingdon OX14 3NB, UK
TEL: +44 (0)1235 529449 FAX: +44 (0)1235 533420
E-MAIL: info@RnDSystems.co.uk

China | R&D Systems China Co., Ltd.

24A1 Hua Min Empire Plaza, 726 West Yan An Road, Shanghai PRC 200050
TEL: +86 (21) 52380373 FAX: +86 (21) 52371001
E-MAIL: info@RnDSystemsChina.com.cn

MATERIALS PROVIDED & STORAGE CONDITIONS

Store the unopened kit at 2-8 °C. Do not use past kit expiration date.

COMPONENTS	PART #	# VIALS	STORAGE OF OPENED/ RECONSTITUTED MATERIAL
Mouse Anti-Human CD3	967554	1 vial	May be stored at 2-8 °C under sterile conditions for up to 30 days or at -20 °C to -70 °C in a manual defrost freezer for up to 3 months.*
Human Th1 Reagent 1	967555	1 vial	
Human Th1 Reagent 2	967556	1 vial	
Reconstitution Buffer 1	967552	1 vial	May be stored under sterile conditions for up to 3 months at 2-8 °C.*
Reconstitution Buffer 2	967553	1 vial	
20X Wash Buffer	967557	3 vials	

* Provided this is within the expiration date of the kit.

OTHER MATERIALS & SUPPLIES REQUIRED

- MagCelect™ Human Naïve CD4⁺ T Cell Isolation Kit (R&D Systems, Catalog # MAGH115, or equivalent).
- Monesin
- PMA
- Ficoll-Hypaque™
- RPMI 1640
- L-Glutamine
- Ionomycin
- Penicillin
- Streptomycin
- Fetal Bovine Serum (FBS)
- β-Mercaptoethanol (2-ME)
- Tissue culture flasks and/or plates
- Microscope
- Hemocytometer
- 37 °C, 5% CO₂ incubator
- Centrifuge
- Pipettes and pipette tips

REAGENT PREPARATION

Human Th1 Differentiation Media

1. Reconstitute Human Th1 Reagent 1 with 250 µL of Reconstitution Buffer 1, this is a 200X stock.
2. Reconstitute Human Th1 Reagent 2 with 250 µL of Reconstitution Buffer 1, this is a 200X stock.
3. Add 50 µL of Human Th1 Reagent 1 and 50 µL of Human Th1 Reagent 2 to 9.9 mL of cell culture media (RPMI, 2 mM L-Glutamine, 50 units/mL Penicillin, 50 µg/mL Streptomycin, 5% FBS, and 50 µM 2-ME).

Human CD3 Antibody

1. Reconstitute the Mouse Anti-Human CD3 antibody with 150 µL of Reconstitution Buffer 2, this is a 100X stock.
2. Add a 10 mL of 20X Wash Buffer to 190 mL of sterile deionized water to prepare 200 mL of 1X Wash Buffer.
3. Just before coating, dilute the 100X antibody stock 100-fold with 1X Wash Buffer.

PROTOCOL FOR Th1 DIFFERENTIATION

1. Coat a plate with Mouse Anti-Human CD3 antibody.
 - a. For a 24-well plate, add 250 µL/well of diluted CD3 antibody. For a 96-well plate, add 50 µL/well of diluted CD3 antibody.
 - b. Incubate at 2-8 °C overnight.
 - c. Wash the plate with 1X Wash Buffer twice before use.
2. Isolate human peripheral blood mononuclear cells (PBMCs) from human blood using Ficoll-Hypaque density gradient centrifugation.
3. Isolate human naïve CD4⁺ T cells from human PBMCs using the MagCelect Human Naïve CD4⁺ T Cell Isolation Kit.
4. Suspend human naïve CD4⁺ T cells at 1-2 x 10⁵ cells/mL in Human Th1 Differentiation Media.
5. Add the cells to a human CD3 antibody-coated plate. For a 24-well plate, add 1 mL/well. For a 96-well plate, add 0.2 mL/well.
6. Incubate the cells in a 37 °C, 5% CO₂ humidified incubator for 5 days.
7. Collect media to analyze cytokine production profile.
8. Wash the cells once with RPMI, resuspend the cells in 1 mL of RPMI, 2 mM L-Glutamine, 50 units/mL penicillin, 50 µg/mL streptomycin, 10% FBS, 50 ng/mL PMA, and 1 µg/mL ionomycin. Incubate the cells in a 37 °C, 5% CO₂, humidified incubator for 1 hour. Then add monesin at 3 µM and incubate for 3 hours.
9. Analyze cytokine expression via flow cytometry.

PROTOCOL OUTLINE

Coat wells of a plate with Mouse Anti-Human CD3 Antibody.

Isolate PBMCs from human blood.

Isolate human naïve CD4⁺ T cells from PBMCs (e.g., using magnetic cell selection).

Perform a cell count.

Suspend 1-2 x 10⁵ naïve CD4⁺ T cells/mL in differentiation media supplemented with Human Th1 Reagents 1 and 2. **Culture** the cells on plates pre-coated with CD3 antibody.

Verify Th1 cell differentiation by analyzing cytokine expression via flow cytometry.

