DIASTAT[®]

Anti-Thyroglobulin (Tg)

For professional use only
Usage reserve aux professionnels
Sólo para uso profesional
Nur für den fachgebrauch
Solo per uso professionale
Endast för professionell användning

Document No. E-23-0106-08
January, 2014

DIASTAT[®] Anti-Thyroglobulin (Tg)

English:	page	2
Français:	page	14
Español:	página	27
Deutsch:	Seite	40
Italiano:	pagina	53
Svenska:	sida	65

REF

FATG 200

IVD

ENGLISH: INTENDED USE

The DIASTAT[®] anti-thyroglobulin (anti-Tg) test is a quantitative/qualitative enzyme-linked immunosorbent assay (ELISA) for the detection of autoantibodies specific for thyroglobulin in human serum or EDTA, Li-Heparin or citrated plasma. It is intended to aid in the diagnosis of autoimmune thyroid disorders and is not definitive in isolation. Autoantibody levels represent one parameter in a multicriterion diagnostic process.

INTRODUCTION

Thyroglobulin (Tg), thyroid peroxidase (TPO) and the TSH receptor are considered the major autoantigens associated with chronic autoimmune thyroiditis. Clinical diagnosis is usually based on the presence of serum or plasma autoantibodies to Tg and TPO¹. In Graves' disease (chronic primary hyperthyroidism), the TSH receptor is the antigen most directly involved in clinical manifestations of the disease².

Thyroglobulin, a 660kD water-soluble glycoprotein³, is the precursor of the thyroid hormones tri-iodothyronine and thyroxine. The presence of anti-Tg autoantibodies in patients with Hashimoto's thyroiditis was first demonstrated in 1956 by Roitt et al⁴ using gel diffusion precipitation. Anti-Tg autoantibodies are predominantly of the IgG class⁵. They are detected, often with TPO autoantibodies, in the majority of Graves' disease and Hashimoto's, including other variants of chronic primary hypothyroidism such as myxoedema and asymptomatic thyroiditis. They have also been detected in spontaneous or post-partum painless thyroiditis⁶, in thyroid autoimmunity with rheumatoid arthritis^{7,8}, and in non-thyroid autoimmune diseases such as Addison's disease and Type 1 diabetes mellitus⁸. Depending upon methodology and assay cut-off, they are also detected in 1.3% to 14.6%^{8,9} of apparently healthy asymptomatic subjects. The role of circulating Tg autoantibodies is unclear; they may simply be indicators of disease as, unlike TPO autoantibodies, they do not appear to be pathogenic.

PRINCIPLE OF THE ASSAY

The wells of the microtitre strips are coated with purified human Tg antigen. During the first incubation, specific autoantibodies in diluted serum or plasma bind to the antigen-coated surface. The wells are then washed to remove unbound components. In the second incubation, the Conjugate, enzyme-labelled antibodies to human IgG, binds any surface-bound autoantibodies. After further washing, specific autoantibodies are traced by incubation with the Substrate. Addition of Stop Solution terminates the reaction, resulting in a coloured end-product. The amount of Conjugate bound is measured in absorbance units. In the qualitative protocol, the amount of Conjugate bound by the sample is compared with that bound by the Reference Control. In the quantitative protocol, the concentration of anti-Tg autoantibody can be estimated by interpolation from a dose-response curve based on Standards. The Standards are referenced against the International Reference Preparation of anti-thyroglobulin 65/93.

KIT COMPONENTS

A	IgG Conjugate 1	1 x 15 mL	Alkaline phosphatase-labelled antibodies to human IgG, Tris buffer, protein stabiliser, <0.1% (w/v) sodium azide. Ready-to-use.	
B	Substrate	1 x 15 mL	Mg ²⁺ , phenolphthalein monophosphate (PMP), buffer solution. Ready-to-use. Do not expose to light during storage.	
C	Stop Solution	1 x 15 mL	Sodium hydroxide, EDTA, carbonate buffer (pH >10). Ready-to-use.	
D	Wash Buffer Concentrate (16X)	3 x 25 mL	Borate buffer, 0.4% (w/v) sodium azide. Dilute before use.	
E	Tg-Coated Wells and Strip Holder	12 x 8 well microtitre strips	Coated with Tg antigen, in a resealable foil pack with desiccant. Colour-coded ROYAL BLUE . Individual wells can be broken off from each microtitre strip.	
F	Sample Diluent Concentrate (5X)	1 x 25 mL	Phosphate buffer, protein stabiliser, 0.5% (w/v) sodium azide. Dilute before use.	
1-5	Anti-Tg Standards	5 x 1.0 mL	Human plasma, buffer, <0.1% (w/v) sodium azide. 0, 100, 400, 1500, 5000 IU/mL. Ready-to-use.	
6	Anti-Tg Reference Control	1 x 1.5 mL	Human plasma, buffer, <0.1% (w/v) sodium azide. Ready-to-use.	
+/-	Positive Control Negative Control	1 x 0.2 mL 1 x 0.1 mL	Human plasma, <0.1% (w/v) sodium azide. Dilute 1:101 with diluted Sample Diluent before use, as for samples.	
	Pack Leaflet			

STORAGE OF REAGENTS

Opened Kit Stability

A kit was opened and reused on three occasions over a three month period with no adverse effect on kit performance.

Handling and Procedural Notes

1. Store kit components at 2-8° C and use until the expiry date on the labels. Do not use expired reagents.
2. Do not mix different lot numbers.
3. Do not freeze kits.
4. Wash Buffer Concentrate, Sample Diluent Concentrate and Positive and Negative Controls must be diluted before use. All other reagents are ready-to-use.
5. Diluted Wash Buffer and diluted Sample Diluent are stable at 2-8° C for up to 6 months if microbial contamination is avoided.
6. Replace surplus microtitre strips in the foil pack and store with the desiccant at 2-8° C, until required.
7. The plate holder is adapted for use with snappable wells only.
8. Do not expose Substrate to light during storage.
9. Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.

Indications of Deterioration

The Substrate should be pale yellow in colour. Pink colouring indicates contamination and the reagent must be discarded. Turbidity or precipitation in any component indicates deterioration and the component should be discarded.

Sample Collection and Storage

The assay is recommended for serum/EDTA, Li-Heparin or citrated plasma samples; do not use lipaemic, haemolysed or turbid samples. Thoroughly mix thawed samples before assay and avoid repeated freeze/thawing. Do not heat-inactivate samples, this may yield false positive results. Samples may be stored undiluted or at 1:101 dilution in diluted Sample Diluent at -20° C or 2-8° C for two weeks.

WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use only.

Safety Precautions

1. Adhere strictly to the instructions in this booklet, particularly for handling and storage conditions.
2. Standards and Controls contain human plasma tested by FDA-cleared assays for hepatitis B surface antigen, HCV, HIV antigen and HIV antibodies and found to be non-reactive/negative. As no known test offers complete assurance that infectious agents are absent, Standards and Controls should be considered potentially infectious and handled with the same precautions as any other potentially biohazardous material. The CDC/NIH Health Manual "Biosafety in Microbiological and Biomedical Laboratories", 3rd edition, 1993, describes how these materials should be handled in accordance with Good Laboratory Practice. This is applicable in the USA.
3. Do not pipette by mouth.
4. Do not smoke, eat, drink or apply cosmetics in areas where kits and samples are handled.
5. Any skin complaints, cuts, abrasions and other skin lesions should be suitably protected.
6. The Standards, Controls, Conjugate, Sample Diluent Concentrate and Wash Buffer Concentrate contain sodium azide which can react with lead and copper plumbing to form highly explosive metal azides. On disposal, drain with large quantities of water to prevent azide build-up.
7. The Stop Solution contains sodium hydroxide. Avoid contact with skin, eyes and mucous membranes. Spillage should be mopped up with copious amounts of water. If contact with skin or eyes occurs, irrigate with water and seek medical attention immediately.
8. The substrate contains PMP, Bronidox L and Diethanolamine. Avoid contact with skin, eyes and respiratory system. If contact with skin, eyes or respiratory system occurs, rinse with water and seek medical advice.
9. Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

Warning

B.

SUBS

Contains: Diethanolamine

H319:	Causes serious eye irritation.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313:	If eye irritation persists: Get medical advice/attention.

C.

SOLN	STOP
------	------

Warning

Contains: Sodium hydroxide

- H315: Causes skin irritation.
 H319: Causes serious eye irritation.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P302+P352: IF ON SKIN: Wash with plenty of soap and water.
 P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P332+P313: If skin irritation occurs: Get medical advice/attention.
 P337+P313: If eye irritation persists: Get medical advice/attention.

D. and F.

BUF	WASH	16X
DIL	SPE	5X

Warning

Contains: Sodium azide

- H302: Harmful if swallowed.
 EUH032: Contact with acids liberates very toxic gas.
 H412: Harmful to aquatic life with long lasting effects.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
 P273: Avoid release to the environment.

PREPARATION

Materials/Equipment Required but not Provided

1. 96 well plate/strip reader with 550nm filter (540-565nm is acceptable).
2. Precision pipettes to dispense 10µL, 100µL, 1mL. Automatic pipette to dispense 100µL. Automatic pipette to dispense 200µL for manual washing, automatic plate washer optional.
3. Glass/plastic measuring cylinders: 1x100mL, 1x400mL.
4. 1mL volume containers.
5. Distilled/deionised water.
6. Paper towels.
7. Timer for 30 and 60 minute intervals.

Preparation for the Assay

Allow all kit components, including the microtitre strips, to warm up to 18-25° C for 30-60 minutes before use. Mix reagents by gentle inversion.

Do not dilute the Reference Control.

Dilute the following reagents and mix thoroughly.

Reagent	Volume	Add
Wash Buffer Concentrate	1 vial	375 mL distilled/deionised water
Sample Diluent Concentrate	1 vial	100 mL distilled/deionised water
Positive and Negative Controls/samples	10µL	1 mL diluted Sample Diluent

Microtitre wells are supplied in strips of eight. If other than a multiple of eight wells are required, proceed as follows.

1. Remove strip from holder by pushing underside of wells.
2. Snap off required number of wells.
3. Hinge rectangular hole into bottom edge (to H) of the holder groove.
4. Ensure the square hole, with nick on left, is firmly held along the top edge (row A).

ASSAY PROTOCOL

Qualitative protocol: run Reference Control, Positive and Negative Controls, and samples.

Quantitative protocol: run Standards (1-5), Positive and Negative Controls, and samples.

1. Reference wells for identification.
2. Pipette 100µL Reference Control/Standards in duplicate, pre-diluted Positive and Negative Controls, and pre-diluted patient samples into appropriate wells. Remember to change pipette tips between additions. This step should not exceed 15 minutes for any one set of Standards/Controls/samples.
3. Incubate 60±10 minutes at 18-25° C.
4. Decant strip contents by quick inversion over a sink suitable for the disposal of biological materials, bearing in mind the potential infective hazard of the samples. Blot inverted strips well with paper towels.
5. Wash wells **five times** with a minimum of 200µL diluted Wash Buffer. **Decant and blot after each wash step.**
6. Add 100µL IgG Conjugate 1 to each well.
7. Incubate 30±5 minutes at 18-25° C.
8. Repeat steps 4 and 5.
9. Add 100µL Substrate to each well.
10. Incubate 30±5 minutes at 18-25°C. **Do not decant.**
11. Add 100µL Stop Solution to each well, in the same order and rate as the Substrate. Tap wells gently to mix.
12. Read strips within 24 hours at 550nm (540-565nm).

CALCULATION AND INTERPRETATION OF RESULTS

Consider each assay separately when calculating and interpreting results.

Qualitative Protocol

Calculate the absorbance value (optical density) ratio for the Positive and Negative Controls, and for each sample.

$$\text{Absorbance Ratio} = \frac{\text{Sample or Control Absorbance Value}}{\text{mean Reference Control Absorbance Value}}$$

Users should calculate a cut-off between positive and negative samples that is specific to their patient populations. Results from the patient populations used in the Euro Diagnostica clinical trial suggest the following cut-off:

Absorbance Ratio

<0.95

≥0.95 to ≤1.0

>1.0

Result Interpretation

Negative

Borderline - recommend repeat testing

Positive

Quantitative Protocol

Plot the mean absorbance value of each Standard against \log_{10} Standard concentration (see following table) on suitable graph paper. Concentrations of Controls and samples can then be read from the standard curve; a typical plot is shown below for reference purposes, it must not be used for interpreting results. 4-parameter logistic (4PL), log/logit or spline curve fits are also satisfactory.


Samples with absorbances above Standard 5 (5000 IU/mL) are outside the range of the assay, and should be reported as >5000 IU/mL, diluted and re-assayed, correcting for this further dilution factor.

NB: As in any assay measuring antibodies, this assay determines the activity of the antibody present in the sample, rather than the concentration. Activity can be affected by a number of parameters, such as antibody avidity.

Standard Concentrations

Standard Number	Concentration IU/mL
1	0
2	100
3	400
4	1500
5	5000

Typical Standard Curve

QUALITY CONTROL

Ensure that adequate maintenance and calibration of the plate-reader is performed according to the manufacturer's instructions, and that the correct wavelength is employed.

Users should ensure they are fully acquainted with the instructions for the assay, particularly the Warnings and Precautions section, and the Handling and Procedural Notes. Users should demonstrate that they can obtain performance specifications for precision and reportable range of test results comparable to those established by the manufacturer before reporting patient test results. It is recommended that the pre-diluted Positive and Negative Controls are run in duplicate in all assays to monitor the quality of the test procedure. Run the ready-to-use Reference Control in duplicate in all qualitative assays.

Assuming the precision specifications described by the manufacturer are met, failure of any Control to meet the Control ratio specifications below renders the assay invalid and patient results should not be reported. The operator may repeat the assay, having reviewed their procedure, or contact the distributor/manufacturer. If repeating the assay, prepare a fresh dilution of each Control and sample. Laboratories may wish to include in-house controls in each assay run. Store such control material at or below -20° C and avoid repeat freeze/thaw cycles. Preservatives such as sodium azide at 0.1% (w/v) will not affect sample results.

Levels of analytes identified in particular diseases are those established by the manufacturer for specific populations, and may not necessarily mirror the literature.

Incidence levels, their relationship to specific diseases, reference ranges, and appropriate cut-off points should all be calculated for the specific populations serviced by users.

Control Ratio Specifications

Protocol	Specifications
Qualitative (ratios)	$\frac{\text{Positive Control Absorbance}}{\text{Reference Control Absorbance}}$ see Positive Control label
	$\frac{\text{Negative Control Absorbance}}{\text{Reference Control Absorbance}} < 0.95$
Quantitative	See Positive Control label for acceptable expected range (IU/mL)
	Negative Control concentration <100 IU/mL

EXPECTED VALUES

207 serum samples from asymptomatic apparently healthy donors, comprising approximately equal numbers of males and females aged between 17 and 68 years with no history of autoimmune disease, were assayed for anti-Tg antibodies of the IgG class. Three samples gave results of 582, 997, and 1396 IU/mL, these were excluded from the cut-off calculations. Using the mean value of 204 samples plus 2.5 standard deviations, the reference range was established with a cut-off of 100 IU/mL. 195/204 (95.6%) of this population gave values less than 100 IU/mL. The distribution of the 4.4% which gave values greater than 100 IU/mL is given below.

<i>Reference Range</i>
≤100 IU/mL = Negative
>100 IU/mL = Positive

No.	Anti-Tg IU/mL	% Asymptomatic Population (n=204)
6	100-119.9	2.9
-	120-139.9	-
2	140-159.9	0.97
1	160-179.9	0.48

Using this reference range, the autoantibody distribution in 116 samples from thyroid disease states was obtained. The diagnostic information is based on the putative diagnosis given by the laboratories that supplied the samples. Patients who tested positive for Paul Bunnel viraemia are also included. The distribution of results is given in the following table.

Putative Diagnosis	n	Range (IU/mL)			
		0 – 100	>100–1000	>1000–5000	>5000
Graves/Thyrotoxicosis	59	35	19	5	-
Hypothyroidism	24	11	10	3	-
Hashimoto's/Thyroiditis/Autoimmune Thyroid Disease	33	25	8	-	-
Paul Bunnel Viraemia Positive	31	30	1	-	-

PERFORMANCE DATA

Agreement with Agglutination Anti-Tg Assay

A comparative study was made with a commercially available anti-Tg agglutination test using 299 patient samples. 28 samples gave a borderline result in the agglutination assay; these were excluded from the overall agreement. The following results were obtained:

		DIASTAT	
		+ve	-ve
AGGLUTINATION ANTI-Tg	+ve	53	5*
	-ve	12≅	201

Overall agreement = 93.7%. The 5* and 12≅ discrepant samples were tested in two commercially available anti-Tg ELISAs. The 5* were negative in both of these ELISAs; of the 12≅, ten were positive in one ELISA and four were positive in the second ELISA.

Agreement with ELISA Anti-Tg Assay

A comparative study was made with another anti-Tg ELISA using 169 patient samples. The following results were obtained:

		DIASTAT	
		+ve	-ve
ELISA ANTI-Tg	+ve	79	6
	-ve	1	83

Overall agreement = 95.9%. The seven discrepant samples were tested in a further ELISA kit. The DIASTAT result was confirmed in 6/7 (85.7%).

Dilution Characteristics

Four dilutions of two patient samples were assayed using two kit batches. The following table shows the mean values obtained and the dilution-corrected recovery.

Sample	Dilution	Mean Value IU/mL	Dilution Corrected % Recovery
1	A	1564	100
	A/2	774	99
	A/4	416	106
	A/8	232	119
2	A	1265	100
	A/2	681	108
	A/4	358	113
	A/8	199	126

Imprecision

- Intra-assay imprecision** determined by testing four controls in 12 assays, using three operators and three kit batches, with replication of four.

Control	Mean Value IU/mL	Root Mean Square %CV
1	195	5.3
2	609	7.8
3	752	6.4
4	1755	10.5

- Inter-assay imprecision** determined by testing four controls in 12 assays, using three operators and three kit batches, with replication of four.

Control	Mean Value IU/mL	SD	%CV
1	195	20.7	10.6
2	609	54.2	8.9
3	752	59.4	7.9
4	1755	254.5	14.5

LIMITATIONS OF USE

1. Although the presence of high titres of antibodies to thyroglobulin is indicative of autoimmune thyroid disease, the data must be considered in light of other clinical and laboratory findings.
2. Some individuals may have high levels of anti-Tg antibodies with little or no evidence of clinical disease. By contrast, some patients with autoimmune thyroid disease may have undetectable levels of these antibodies.
3. Low anti-Tg titres may be found in apparently healthy individuals. The clinical significance of this information is currently unclear.
4. For repeat patient sampling, e.g. for monitoring, the same type of sample (serum or EDTA, Li-Heparin or citrated plasma) should be used throughout the study period.

REFERENCES

1. Kohno Y, et al. *J Clin Endocr Met*, **67** No. 5, 899-907, 1988.
2. Rapaport B. *Exp Clin Endocrinol*, **97** No. 2/3, 147-152, 1991.
3. Dorrein Y, et al. *Biochem Biophys Acta*, **2**, 454, 1948.
4. Roitt IM, et al. *Lancet II*, 820-821, 1956.
5. Delespesse G, et al. *Horm Metabl Res*, **8**, 50-54, 1976.
6. Amino N, et al. *J Clin Endoc Metab*, **53**, 113-116, 1981.
7. Walker DJ, et al. *Ann Rheum Dis*, **45**, 323-326, 1986.
8. Scherbaum WA, et al. *Acta Endocrinol*, (Copenh) Supple. 281, 325-329, 1987.
9. Shimijo N, et al. *Clinica Acta*, **163**, 41-49, 1987.

SUMMARY OF PROTOCOL

1. Dilute samples and Positive and Negative Controls 1:101. Do not dilute Standards or Reference Control.
2. Add 100µL of Reference Control/Standards in duplicate, pre-diluted Positive and Negative Controls and samples into referenced wells of the microtitre strip.
3. Incubate 60±10 minutes at 18-25° C.
4. Wash strips 5 times.
5. Add 100µL of IgG Conjugate 1 to each well.
6. Incubate 30±5 minutes at 18-25° C.
7. Wash strips 5 times.
8. Add 100µL of Substrate to each well.
9. Incubate 30±5 minutes at 18-25° C.
10. Add 100µL of Stop Solution to each well.
11. Read absorbance at 550nm.

	Batch code / Code du lot / Código de lote / Chargenbezeichnung / Codice del lotto / Satsnummer
	Catalogue number/ Référence du catalogue / Número de catálogo / Bestellnummer / Numero di catalogo / Katalognummer
	Use by/ Utiliser jusque / Fecha de caducidad / Verwendbar bis / Utilizzare entro / Hållbar till
	Temperature limitation/ Limites de temperature / Limite de temperatura / Zulässiger Temperatur-bereich / Limite di temperatura / Temperaturgränser
	Biological risks/ Risques biologiques / Riesgo biológico/Biogefährdung /Rischio biologico / Biologisk risk
	Consult instructions for use/ Consulter les instructions d'utilisation / Consulte las instruccione de uso/ Gebrauchsanweisung beachten / Consultare le istruzioni per l'uso / Se bruksanvisning
	In Vitro Diagnostic Medical Device/ Dispositif medical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In Vitro Diagnostikum / Dispositivo medico-diagnostico in vitro / In vitro diagnostika
	Manufacturer/ Fabricant / Fabricante / Hersteller / Fabbicante / Tillverkare
	Contains sufficient for <n> tests/ Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos/ Ausreichend für "n" Ansätze / Conteúdo suficiente para <n> ensaios / Innehållet tillräckligt för n tester
	Warning / Attention / Atención / Achtung / Attenzione / Varning
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive / Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro / Konform mit Richtlinie 98/79/EG zu In-vitro- Diagnostika / Conformità alla direttiva 98/79/CE relativa ai dispositivi medico- diagnostici in vitro / Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

CONJ	Conjugate/ Conjugué / Conjugado / Konjugat / Conjugato / Konjugat
SUBS	Substrate/ Substrat / Substrato / Substrat / Substrato / Substrat
SOLN STOP	Stop solution/ Solution d'Arrêt / Solución de Parada / Stopplösung/ Soluzione bloccante / Stopplösning
BUF WASH 16 x	Wash buffer concentrate (16 X)/ Concentré tampon de (16X lavage)/ Concentrado de Búfer de lavado (16 X) / Washpuffer-Konzentrat (16 X) / Tampone di lavaggio concentrato (16 X) / Tvättbuffert koncentrat (16 X)
Ag	Tg-coated wells and strip holder / Cupules enduites de Tg et Portebandes / Soporte para Bandas y Vasos Recubiertos con Tg / Tg-beschichtete Vertiefungen und Streifenrahmen / Pozzetti rivestiti di Tg e supporto per strip /Tg-klädda brunnar och striphållare
DIL SPE 5 X	Sample Diluent Concentrate (5 X) / Concentré diluent pour échantillons (5 X)/ Concentrado de Diluente de Muestra (5 X) / Probendiluens Konzentrat / Diluente per campioni concentrato (5 X) / Provspädningsbuffert koncentrat (5 X)
CAL X	Anti-Tg Standards 1-5/ Etalons anti-Tg 1-5 / Estandares Anti-Tg 1-5 / Anti-Tg Standards 1-5 / Standard anti-Tg 1-5 / Anti-Tg-standarder
CONTROL REF	Anti-Tg Reference Control/ Témoins de référence anti-Tg / Control de Referencia Anti-Tg / Anti-Tg Referenzkontrolle / Controllo di riferimento anti-Tg / Anti-Tg referenskontroll
CONTROL +	Positive Controls/ Témoins positives / Controles Positivos / Positiv-Kontrollen / Controlli Positivi / Positiva kontrollen
CONTROL -	Negative Controls/ Témoins négatifs / Controles Negativos / Negativ-Kontrollen / Controlli negativi / Negativa kontrollen

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EURO DIAGNOSTICA AB