DIASTAT[®]

Anti-Mitochondria

For professional use only
Usage reserve aux professionnels
Sólo para uso profesional
Nur für den fachgebrauch
Solo per uso professionale
Endast för professionell användning

Document No. E-23-0117-06
January, 2014

DIASTAT[®] Anti-Mitochondria

English: page 2
Français: page 11
Español: página 21
Deutsch: Seite 31
Italiano: pagina 40
Svenska: sida 49

REF

FMIT 200

IVD

ENGLISH : INTENDED USE

The DIASTAT[®] Anti-mitochondrial antibody test is a qualitative enzyme-linked immunosorbent assay (ELISA) for the detection of IgG and IgM autoantibodies specific for the M2 antigen of the mitochondrial membrane in human serum or plasma. Detection of these autoantibodies is a useful aid in the diagnosis of Primary Biliary Cirrhosis and is not definitive in isolation. Autoantibody levels represent one parameter in a multicriterion diagnostic process.

INTRODUCTION

Primary Biliary Cirrhosis (PBC) is a chronic and often fatal cholestatic autoimmune liver disease, primarily affecting young to middle-aged women^{1,2}. It is characterised by slow, progressive destruction of intrahepatic bile ducts, portal inflammation and scarring, eventually leading to cirrhosis and liver failure. Clinical manifestations, biochemical tests and the detection of circulating anti-mitochondrial antibodies (AMAs) are used for diagnosis of PBC. Histological examination provides confirmatory information and is helpful in assessing prognosis. AMAs are reportedly present in over 95% of biopsy-proven PBCs^{3,4}, they have been observed at low levels in patients with other diseases^{4,5} and in apparently healthy individuals^{6,7}. Disease-specific subtypes of AMAs have been identified and are regarded as specific and sensitive markers for PBC. To date, several autoantigens have been identified⁸; of these, the M2 antigen is believed to be specific for PBC^{4,9,10}. It is composed of two immunoreactive polypeptides of approximately 70kD and 50kD and has been identified as part of the pyruvate dehydrogenase complex, a major mitochondrial enzyme^{4,10}.

Current techniques available for the detection of AMAs include immunofluorescence, radio-immunoassay, enzyme-immunoassay and immunoblotting. The most commonly used method is immunofluorescence on the HEp 2 substrate, however immunofluorescent staining patterns can be non-specific, making identification of the antigen subtype difficult. ELISAs offer increased sensitivity, simplicity and more rapid results.

PRINCIPLE OF THE ASSAY

The wells of the microtitre strips are coated with M2 antigen. During the first incubation, specific autoantibodies in diluted serum or plasma bind to the antigen-coated surface. The wells are then washed to remove unbound components. In the second incubation, the Conjugate, enzyme-labelled antibodies to human IgG and IgM, binds any surface-bound autoantibodies. After further washing, specific autoantibodies are traced by incubation with the Substrate. Addition of Stop Solution terminates the reaction, resulting in a coloured end-product. The amount of Conjugate bound is measured in absorbance units and compared with that bound by the Reference Control.

KIT COMPONENTS

A	IgG/IgM Conjugate	1 x 15mL	Alkaline phosphatase-labelled antibodies to human IgG and IgM, Tris buffer, protein stabiliser, <0.1% (w/v) sodium azide. Ready-to-use.	
B	Substrate	1 x 15mL	Mg ²⁺ , phenolphthalein monophosphate (PMP), buffer solution. Ready-to-use. Do not expose to light during storage.	
C	Stop Solution	1 x 15mL	Sodium hydroxide, EDTA, carbonate buffer (pH >10). Ready-to-use.	
D	Wash Buffer Concentrate (16X)	2 x 25mL	Borate buffer, 0.4% (w/v) sodium azide. Dilute before use.	
E	M2-Coated Wells and Strip Holder	12 x 8 well microtitre strips	Coated with M2 antigen, in a resealable foil pack with desiccant. Colour-coded PINK . Individual wells can be broken off from each microtitre strip.	
F	Sample Diluent Concentrate (5X)	1 x 25mL	Phosphate buffer, protein stabiliser, 0.5% (w/v) sodium azide. Dilute before use.	
6	Anti-M2 Reference Control	1 x 1.5mL	Human plasma, buffer, <0.1% (w/v) sodium azide. Ready-to-use.	
+/-	Positive Control Negative Control	1 x 0.2mL 1 x 0.1mL	Human plasma, <0.1% (w/v) sodium azide. Dilute 1:101 with diluted Sample Diluent before use, as for samples.	
	Pack Leaflet			

STORAGE OF REAGENTS

Opened Kit Stability

A kit was opened and reused on three occasions over a three month period with no adverse effect on kit performance.

Handling and Procedural Notes

1. Store kit components at 2-8° C and use until the expiry date on the labels. Do not use expired reagents.
2. Do not mix different lot numbers.
3. Do not freeze kits.
4. Wash Buffer Concentrate, Sample Diluent Concentrate and Positive and Negative Controls must be diluted before use. All other reagents are ready-to-use.
5. Diluted Wash Buffer and diluted Sample Diluent are stable at 2-8° C for up to 6 months if microbial contamination is avoided.
6. Replace surplus microtitre strips in the foil pack and store with the desiccant at 2-8° C, until required.
7. The plate holder is adapted for use with snappable wells only.
8. Do not expose Substrate to light during storage.
9. Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.

Indications of Deterioration

The Substrate should be pale yellow in colour. Pink colouring indicates contamination and the reagent must be discarded. Turbidity or precipitation in any component indicates deterioration and the component should be discarded.

Sample Collection and Storage

The assay is recommended for serum/plasma samples; do not use lipaemic, haemolysed or turbid samples. Thoroughly mix thawed samples before assay and avoid repeated freeze/thawing. Do not heat-inactivate samples, this may yield false positive results.

Samples may be stored at or below -20° C, or assayed within 12 hours of collection.

WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use only.

Safety Precautions

1. Adhere strictly to the instructions in this booklet, particularly for handling and storage conditions.
2. The Controls contain human plasma tested by FDA-cleared assays for hepatitis B surface antigen, HCV, HIV antigen and HIV antibodies and found to be non-reactive/negative. As no known test offers complete assurance that infectious agents are absent, Controls should be considered potentially infectious and handled with the same precautions as any other potentially biohazardous material. The CDC/NIH Health Manual "Biosafety in Microbiological and Biomedical Laboratories", 3rd edition, 1993, describes how these materials should be handled in accordance with Good Laboratory Practice. This is applicable in the USA.
3. Do not pipette by mouth.
4. Do not smoke, eat, drink or apply cosmetics in areas where kits and samples are handled.
5. Any skin complaints, cuts, abrasions and other skin lesions should be suitably protected.
6. The Controls, Conjugate, Sample Diluent Concentrate and Wash Buffer Concentrate contain sodium azide which can react with lead and copper plumbing to form highly explosive metal azides. On disposal, drain with large quantities of water to prevent azide build-up.
7. The Stop Solution contains sodium hydroxide. Avoid contact with skin, eyes and mucous membranes. Spillage should be mopped up with copious amounts of water. If contact with skin or eyes occurs, irrigate with water and seek medical attention immediately.
8. The substrate contains PMP, Bronidox L and Diethanolamine. Avoid contact with skin, eyes and respiratory system. If contact with skin, eyes or respiratory system occurs, rinse with water and seek medical advice.
9. Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

Warning

B.

SUBS

Contains: Diethanolamine

- | | |
|-----------------|--|
| H319: | Causes serious eye irritation. |
| P264: | Wash hands thoroughly after handling. |
| P280: | Wear protective gloves/protective clothing/eye protection/face protection. |
| P305+P351+P338: | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313: | If eye irritation persists: Get medical advice/attention. |

C.

SOLN	STOP
------	------

Warning

Contains: Sodium hydroxide

- H315: Causes skin irritation.
 H319: Causes serious eye irritation.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P302+P352: IF ON SKIN: Wash with plenty of soap and water.
 P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P332+P313: If skin irritation occurs: Get medical advice/attention.
 P337+P313: If eye irritation persists: Get medical advice/attention.

D. and F.

BUF	WASH	16X
DIL	SPE	5X

Warning

Contains: Sodium azide

- H302: Harmful if swallowed.
 EUH032: Contact with acids liberates very toxic gas.
 H412: Harmful to aquatic life with long lasting effects.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
 P273: Avoid release to the environment.

P R E P A R A T I O N

Materials/Equipment Required but not Provided

1. 96 well plate/strip reader with 550nm filter (540-565nm is acceptable).
2. Precision pipettes to dispense 10µL, 100µL, 1mL. Automatic pipette to dispense 100µL. Automatic pipette to dispense 200µL for manual washing, automatic plate washer optional.
3. Glass/plastic measuring cylinders: 1x100mL, 1x400mL.
4. 1mL volume containers.
5. Distilled/deionised water.
6. Paper towels.
7. Timer for 30 and 60 minute intervals.

Preparation for the Assay

Allow all kit components, including the microtitre strips, to warm up to 18-25° C for 30-60 minutes before use. Mix reagents by gentle inversion.

Do not dilute the Reference Control.

Dilute the following reagents and mix thoroughly.

Reagent	Volume	Add
Wash Buffer Concentrate	1 vial	375mL distilled/deionised water
Sample Diluent Concentrate	1 vial	100mL distilled/deionised water
Positive and Negative Controls/samples	10µL	1mL diluted Sample Diluent

Microtitre wells are supplied in strips of eight. If other than a multiple of eight wells are required, proceed as follows.

1. Remove strip from holder by pushing underside of wells.
2. Snap off required number of wells.
3. Hinge rectangular hole into bottom edge (to H) of the holder groove.
4. Ensure the square hole, with nick on left, is firmly held along the top edge (row A).

ASSAY PROTOCOL

1. Reference wells for identification.
2. Pipette 100µL Reference Control in duplicate, pre-diluted Positive and Negative Controls, and pre-diluted patient samples into appropriate wells. Remember to change pipette tips between additions. This step should not exceed 15 minutes for any one set of Controls/samples.
3. Incubate 60±10 minutes at 18-25° C.
4. Decant strip contents by quick inversion over a sink suitable for the disposal of biological materials, bearing in mind the potential infective hazard of the samples. Blot inverted strips well with paper towels.
5. Wash wells **three times** with a minimum of 200µL diluted Wash Buffer. **Decant and blot after each wash step.**
6. Add 100µL IgG/IgM Conjugate to each well.
7. Incubate 30±5 minutes at 18-25° C.
8. Repeat steps 4 and 5.
9. Add 100µL Substrate to each well.
10. Incubate 30±5 minutes at 18-25° C. **Do not decant.**
11. Add 100µL Stop Solution to each well, in the same order and rate as the Substrate. Tap wells gently to mix.
12. Read strips within 24 hours at 550nm (540-565nm).

CALCULATION AND INTERPRETATION OF RESULTS

Consider each assay separately when calculating and interpreting results.

Calculate the absorbance value (optical density) ratio for the Positive and Negative Controls, and for each sample.

$$\text{Absorbance Ratio} = \frac{\text{Sample or Control Absorbance Value}}{\text{mean Reference Control Absorbance Value}}$$

Users should calculate a cut-off between positive and negative samples that is specific to their patient populations. Results from the patient populations used in the Euro Diagnostica clinical trial suggest the following cut-off:

Absorbance Ratio
<0.95
≥0.95 to ≤1.0
>1.0

Result Interpretation
Negative
Borderline - recommend repeat testing
Positive

QUALITY CONTROL

Ensure that adequate maintenance and calibration of the plate-reader is performed according to the manufacturer's instructions, and that the correct wavelength is employed.

Users should ensure they are fully acquainted with the instructions for the assay, particularly the Warnings and Precautions section, and the Handling and Procedural Notes. Users should demonstrate that they can obtain performance specifications for precision and reportable range of test results comparable to those established by the manufacturer before reporting patient test results. It is recommended that the pre-diluted Positive and Negative Controls are run in duplicate in all assays to monitor the quality of the test procedure. Run the ready-to-use Reference Control in duplicate in all assays.

Assuming the precision specifications described by the manufacturer are met, failure of any Control to meet the Control ratio specifications below renders the assay invalid and patient results should not be reported. The operator may repeat the assay, having reviewed their procedure, or contact the distributor/manufacturer. If repeating the assay, prepare a fresh dilution of each Control and sample. Laboratories may wish to include in-house controls in each assay run. Store such control material at or below -20° C and avoid repeat freeze/thaw cycles. Preservatives such as sodium azide at 0.1% (w/v) will not affect sample results.

Levels of analytes identified in particular diseases are those established by the manufacturer for specific populations, and may not necessarily mirror the literature. Incidence levels, their relationship to specific diseases, reference ranges, and appropriate cut-off points should all be calculated for the specific populations serviced by users.

Control Ratio Specifications	
Positive Control Absorbance Reference Control Absorbance	see Positive Control label
Negative Control Absorbance Reference Control Absorbance	<1.0

EXPECTED VALUES

Serum samples were collected from reference centres and from asymptomatic healthy individuals. Using the sample to Reference Control ratio of <1 as a cut-off, greater than 99% of the asymptomatic population were negative for anti-mitochondrial IgG and IgM antibodies. The distribution of results found in specific disease groups is given in the following table.

Control and Disease Groups	n	% Positive using Sample to Reference Control ratio > 1
Asymptomatic Healthy	193	0.5
Primary Biliary Cirrhosis	72	95.8
Anti-Centromere Antibody Positive	17	11.8
Systemic Lupus Erythematosus	62	3.2
Chronic Active Hepatitis	5	-
Scleroderma	18	-
Sjögren's Syndrome	21	-
Viral Infections	23	-
Rheumatoid Arthritis	19	-
Polymyositis	23	-

LIMITATIONS OF USE

1. Although the presence of high titres of antibodies to M2 antigen is indicative of Primary Biliary Cirrhosis, the data must be considered in light of other clinical and laboratory findings.
2. Some individuals may have high levels of anti-M2 antibodies with little or no evidence of clinical disease. By contrast, a small number of patients with PBC may have undetectable levels of these antibodies.
3. For repeat patient sampling, e.g. for monitoring, the same type of sample (serum or plasma) should be used throughout the study period.

REFERENCES

1. Kaplan MM. *New England Journal of Medicine*, **316**, 9, 521-529, 1987.
2. Powell FC, Schroeter AL, Dickson ER. *Quarterly Journal of Medicine*, New Series 62, **237**, 75-82, 1987.
3. Kyriatsoulis A, Manns M, et al. *J Immun Methods*, **109**, 113-121, 1988.
4. McHugh NJ, James IE, et al. *Clin Exp Immunol*, **81**, 244-249, 1990.
5. Neuberger J, Lombard M, Galbraith R, *Gut Supplement*, S73-S78, 1991.
6. Berg PA, Homberg JC, et al. *Lancet*, October 10, 1981.
7. Manns M. *Journal of Hepatology Review*, **9**, 272-280, 1989.
8. Surh CD, Danner DJ, et al. *Hepatology*, **9**, 1, 63-68, 1989.
9. Gupta RC. *Clin Exp Immunol*, **59**, 604-612, 1985.
10. Linderborn-Fotinos J, Baum H, Berg PA. *Hepatology*, **5**, 5, 763-769, 1985.

SUMMARY OF PROTOCOL

1. Dilute samples and Positive and Negative Controls 1:101. Do not dilute Reference Control.
2. Add 100µL of Reference Control in duplicate, pre-diluted Positive and Negative Controls and samples into referenced wells of the microtitre strip.
3. Incubate 60±10 minutes at 18-25° C.
4. Wash strips 3 times.
5. Add 100µL of IgG/IgM Conjugate to each well.
6. Incubate 30±5 minutes at 18-25° C.
7. Wash strips 3 times.
8. Add 100µL of Substrate to each well.
9. Incubate 30±5 minutes at 18-25° C.
10. Add 100µL of Stop Solution to each well.
11. Read absorbance at 550nm.

	Batch code / Code du lot / Código de lote / Chargenbezeichnung / Codice del lotto / Satsnummer
	Catalogue number/ Référence du catalogue / Número de catálogo / Bestellnummer / Numero di catalogo / Katalognummer
	Use by/ Utiliser jusque / Fecha de caducidad / Verwendbar bis / Utilizzare entro / Hållbar till
	Temperature limitation/ Limites de temperature / Limite de temperatura / Zulässiger Temperatur-bereich / Limite di temperatura / Temperaturgränser
	Biological risks/ Risques biologiques / Riesgo biológico/Biogefährdung /Rischio biologico / Biologisk risk
	Consult instructions for use/ Consulter les instructions d'utilisation / Consulte las instruccione de uso/ Gebrauchsanweisung beachten / Consultare le istruzioni per l'uso / Se bruksanvisning
	In Vitro Diagnostic Medical Device/ Dispositif medical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In Vitro Diagnostikum / Dispositivo medico-diagnostico in vitro / In vitro diagnostika
	Manufacturer/ Fabricant / Fabricante / Hersteller / Fabbicante / Tillverkare
	Contains sufficient for <n> tests/ Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos/ Ausreichend für "n" Ansätze / Conteúdo suficiente para <n> ensaios / Innehållet tillräckligt för n tester
	Warning / Attention / Atención / Achtung / Attenzione / Varning
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive / Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro / Konform mit Richtlinie 98/79/EG zu In-vitro- Diagnostika / Conformità alla direttiva 98/79/CE relativa ai dispositivi medico- diagnostici in vitro / Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

CONJ	Conjugate/ Conjugué / Conjugado / Konjugat / Conjugato / Konjugat
SUBS	Substrate/ Substrat / Substrato / Substrat / Substrato / Substrat
SOLN STOP	Stop solution/ Solution d'Arrêt / Solución de Parada / Stopplösung/ Soluzione bloccante / Stopplösning
BUF WASH 16 x	Wash buffer concentrate (16 X)/ Concentré tampon de (16X lavage)/ Concentrado de Búfer de lavado (16 X) / Washpuffer-Konzentrat (16 X) / Tampone di lavaggio concentrato (16 X) / Tvättbuffert koncentrat (16 X)
Ag	M2-coated wells and strip holder / Cupules enduites de M2 et Portebandes / Soporte para Bandas y Vasos Recubiertos con M2 / M2-beschichtete Vertiefungen und Streifenrahmen / Pozzetti rivestiti di M2 e supporto per strip /M2-klädda brunnar och striphållare
DIL SPE 5 X	Sample Diluent Concentrate (5 X) / Concentré diluent pour échantillons (5 X)/ Concentrado de Diluente de Muestra (5 X) / Probendiluens Konzentrat / Diluente per campioni concentrato (5 X) / Provspädningsbuffert koncentrat (5 X)
CAL X	Anti-M2 Standards 1-5/ Etalons anti-M2 1-5 / Estandares Anti-M2 1-5 / Anti-M2 Standards 1-5 / Standard anti-M2 1-5 / Anti-M2-standarder
CONTROL REF	Anti-2 Reference Control/ Témoins de référence anti-M2 / Control de Referencia Anti-M2 / Anti-M2 Referenzkontrolle / Controllo di riferimento anti-M2 / Anti-M2 referenskontroll
CONTROL +	Positive Controls/ Témoins positives / Controles Positivos / Positiv-Kontrollen / Controlli Positivi / Positiva kontrollen
CONTROL -	Negative Controls/ Témoins négatifs / Controles Negativos / Negativ-Kontrollen / Controlli negativi / Negativa kontrollen

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EURO DIAGNOSTICA AB