

# **Quick-16S™ Plus NGS Library Prep Kit (V4) - 384 Preps**

Fastest, normalization-free 16S library prep

## Highlights

- **Fast:** 8 times less hands-on time than conventional 16S library prep protocols.
- **Easy:** Premixed Plate minimizes hands-on time and handling errors.
- **Normalization-free:** 100% automation ready with only a single PCR step and no need for normalization. Just pool by equal volume!

Catalog Numbers:  
D6432



Scan with your smart-phone camera to  
view the online protocol/video.



# Table of Contents

---

<b>Product Contents</b> .....	<b>01</b>
<b>Specifications</b> .....	<b>02</b>
<b>Product Description</b> .....	<b>03-04</b>
<b>Protocol</b> .....	<b>05-09</b>
Before Starting .....	<b>05</b>
(I) 1-Step PCR .....	<b>05-06</b>
(II) Pooling by Equal Volume.....	<b>06</b>
(III) Final Library Clean-up .....	<b>07</b>
(IV) Library Quantification .....	<b>08</b>
(V) Denaturation & Dilution.....	<b>09</b>
(VI) Spiking in Sequencing Primers .....	<b>09</b>
<b>Appendices</b> .....	<b>10-11</b>
Removal of PCR Inhibitors from Starting DNA.....	<b>10</b>
Adjusting PCR Cycles .....	<b>10</b>
Sequencing Recommendations for Other Platforms .....	<b>10</b>
Illumina MiSeq Sample Sheet Setup .....	<b>10</b>
Illumina Sequencing Primer Positions .....	<b>11</b>
<b>Complete Your Workflow</b> .....	<b>12</b>
<b>Troubleshooting Guide</b> .....	<b>13</b>
<b>Notes</b> .....	<b>14-16</b>
<b>Guarantee</b> .....	<b>17</b>

# Product Contents

<b>Quick-16S™ Plus NGS Library Prep Kit (V4)</b>	<b>D6432 (384 rxns.)</b>	<b>Storage Temp.</b>
<b>V4 Targeted Amplification Plate<sup>1</sup></b>	8 µl each well	-20°C
<b>Read 1 Sequencing Primer</b>	30 µl	-20°C
<b>Read 2 Sequencing Primer</b>	30 µl	-20°C
<b>Index 1 (i7) Sequencing Primer</b>	30 µl	-20°C
<b>Index 2 (i5) Sequencing Primer<sup>2</sup></b>	30 µl	-20°C
<b>ZymoBIOMICS™ Microbial Community DNA Standard (50 ng)</b>	10 µl	-20°C
<b>Select-a-Size™ MagBead Concentrate<sup>3</sup></b>	30 µl	4-8°C
<b>Select-a-Size™ MagBead Buffer<sup>3</sup></b>	1 ml	4-8°C
<b>DNA Wash Buffer</b>	6 ml	Room Temp.
<b>ZymoBIOMICS™ DNase/RNase Free Water</b>	1 ml	Room Temp.
<b>PCR Inactivation Solution</b>	100 µl	Room Temp.
<b>Magnetic Rods</b>	4	-
<b>Instruction Manual</b>	1 pc	-

<sup>1</sup> Protect the plate from light.

<sup>2</sup> Index 2 primer is included for sequencers that require an i5 index sequencing primer.

<sup>3</sup> The *Select-a-Size*™ MagBead Concentrate and Buffer are shipped at room temperature but should be stored at 4-8°C upon receipt.

# Specifications

- **Sample Input** – Purified microbial DNA ( $\leq 100$  ng), free of PCR inhibitors.
- **V4 Primer Sequences (adapters not included)** – 515f (GTGYCAGCMGCCGCGGTAA) and 806r (GGACTACNVGGGTWTCTAA).
- **Index Sequences** – 10 bp indexes are listed in the MiSeq Sample Sheet Template provided. The template is available for download by visiting the Documentation section of the Product Page at [www.zymoresearch.com](http://www.zymoresearch.com).
- **Amplicon Size** – The final amplicon size after 1-Step PCR (targeted amplification and barcode addition) is  $\sim 388$  bp.
- **Sequencing Platform** – Compatible with all Illumina® sequencing platforms using custom sequencing primers. We recommend the MiSeq® Reagent Kit v2 (300-cycle).
- **Equipment Needed** (user provided) – Microcentrifuge, plate spinner (centrifuge), 384-well real-time quantitative PCR system. The included plate is compatible with the following Applied Biosystems™ PCR systems: ProFlex, Veriti, VeritiPro™, 9700; QuantStudio 5, 6, 6 Pro, 7, 7 Pro, 12K; ViiA 7; 7900 HT.
- **V4 Targeted Amplification Plate** – Each position on the targeted amplification plate contains 5  $\mu$ l of Equalase™ qPCR Premix, 1  $\mu$ l of water, and 1  $\mu$ l of each uniquely indexed (UDI) 515f and 806r V4 primer. Total reaction volume is 10  $\mu$ l upon addition of 2  $\mu$ l DNA as per the protocol.
- **V4 UDI Plate Layout** –

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	Z80001	Z80193	Z80009	Z80201	Z80017	Z80209	Z80025	Z80127	Z80033	Z80025	Z80041	Z80233	Z80049	Z80041	Z80057	Z80049	Z80055	Z80257	Z80073	Z80255	Z80081	Z80271	Z80089	Z80281
B	Z80097	Z80208	Z80055	Z80207	Z80113	Z80095	Z80124	Z80118	Z80129	Z80031	Z80137	Z80039	Z80145	Z80037	Z80153	Z80045	Z80141	Z80053	Z80149	Z80163	Z80177	Z80169	Z80145	Z80177
C	Z80062	Z80194	Z80010	Z80202	Z80018	Z80210	Z80026	Z80018	Z80054	Z80026	Z80042	Z80234	Z80050	Z80042	Z80058	Z80020	Z80066	Z80258	Z80074	Z80266	Z80082	Z80074	Z80090	Z80282
D	Z80098	Z80290	Z80106	Z80298	Z80114	Z80106	Z80122	Z80134	Z80130	Z80022	Z80138	Z80030	Z80146	Z80038	Z80154	Z80046	Z80162	Z80054	Z80170	Z80182	Z80170	Z80390	Z80186	Z80278
E	Z80090	Z80125	Z80013	Z80203	Z80019	Z80111	Z80027	Z80129	Z80099	Z80127	Z80043	Z80129	Z80051	Z80044	Z80099	Z80121	Z80067	Z80259	Z80076	Z80267	Z80083	Z80129	Z80091	Z80219
F	Z80099	Z80291	Z80107	Z80299	Z80115	Z80107	Z80123	Z80135	Z80131	Z80023	Z80139	Z80031	Z80147	Z80039	Z80155	Z80047	Z80163	Z80055	Z80171	Z80183	Z80179	Z80371	Z80187	Z80279
G	Z80004	Z80196	Z80012	Z80024	Z80000	Z80012	Z80028	Z80020	Z80006	Z80028	Z80044	Z80236	Z80052	Z80044	Z80090	Z80022	Z80068	Z80260	Z80076	Z80268	Z80084	Z80276	Z80092	Z80284
H	Z80100	Z80292	Z80109	Z80300	Z80116	Z80106	Z80124	Z80106	Z80112	Z80034	Z80140	Z80032	Z80148	Z80040	Z80156	Z80048	Z80144	Z80056	Z80172	Z80164	Z80166	Z80172	Z80118	Z80039
I	Z80005	Z80197	Z80013	Z80025	Z80021	Z80213	Z80029	Z80021	Z80021	Z80029	Z80045	Z80237	Z80053	Z80045	Z80061	Z80023	Z80069	Z80261	Z80077	Z80269	Z80085	Z80277	Z80093	Z80285
J	Z80011	Z80293	Z80109	Z80301	Z80117	Z80109	Z80125	Z80137	Z80133	Z80025	Z80141	Z80033	Z80149	Z80041	Z80197	Z80049	Z80055	Z80257	Z80173	Z80185	Z80181	Z80373	Z80189	Z80281
K	Z80096	Z80126	Z80014	Z80204	Z80002	Z80124	Z80098	Z80122	Z80098	Z80120	Z80066	Z80138	Z80054	Z80046	Z80062	Z80124	Z80090	Z80262	Z80078	Z80170	Z80068	Z80279	Z80094	Z80286
L	Z80102	Z80294	Z80110	Z80302	Z80118	Z80110	Z80126	Z80118	Z80134	Z80026	Z80142	Z80034	Z80150	Z80042	Z80158	Z80030	Z80066	Z80258	Z80174	Z80166	Z80182	Z80374	Z80190	Z80282
M	Z80087	Z80199	Z80015	Z80207	Z80003	Z80215	Z80031	Z80023	Z80039	Z80021	Z80047	Z80239	Z80055	Z80047	Z80061	Z80025	Z80071	Z80263	Z80079	Z80271	Z80087	Z80279	Z80095	Z80287
N	Z80010	Z80295	Z80113	Z80303	Z80119	Z80111	Z80127	Z80139	Z80115	Z80027	Z80143	Z80035	Z80151	Z80043	Z80159	Z80047	Z80063	Z80259	Z80175	Z80167	Z80183	Z80375	Z80191	Z80283
O	Z80008	Z80000	Z80016	Z80008	Z80024	Z80016	Z80032	Z80024	Z80040	Z80022	Z80048	Z80240	Z80056	Z80048	Z80064	Z80026	Z80072	Z80084	Z80090	Z80272	Z80088	Z80280	Z80096	Z80288
P	Z80104	Z80296	Z80112	Z80304	Z80120	Z80112	Z80128	Z80130	Z80136	Z80028	Z80144	Z80036	Z80152	Z80044	Z80160	Z80052	Z80168	Z80190	Z80176	Z80168	Z80184	Z80376	Z80192	Z80384


Primer Plate A  
 Primer Plate B  
 Primer Plate C  
 Primer Plate D

# Product Description

16S rRNA gene sequencing is a routine technique for microbiome composition profiling. Compared to shotgun metagenomics sequencing, 16S rRNA gene sequencing is more cost-effective and more robust; it generally requires less input DNA and is less impacted by the presence of non-microbial DNA. However, 16S rRNA gene sequencing has its own challenges. Common 16S library preparation protocols have not been optimized to be cost-effective for large-scale applications.


The **Quick-16S™ Plus NGS Library Prep Kit (V4)** is the fastest and simplest library prep method for high-throughput 16S rRNA sequencing. Distinguishing features of the kit are described below.

**Fastest 16S Workflow.** The **Quick-16S™ Plus NGS Library Prep Kit (V4)** utilizes a single qPCR/PCR for combined targeted amplification and barcode addition using specially designed primers. After pooling by equal volume, a single clean-up of the final library is performed, rather than multi-well magnetic bead clean-ups. Additional library quantification analysis such as TapeStation® analysis or gel electrophoresis is not necessary. The addition of custom sequencing primers allows the library to be sequenced on any Illumina 300 cycle kit and ensures there is sufficient overlap between paired reads. With this workflow, the hands-on time of 16S library preparation is reduced to only 30 minutes (Figure 1).




**Figure 1. Quick-16S™ Plus NGS Library Prep Kit (V4) workflow versus the Illumina Metagenomic Protocol.** Total hands-on time calculations are based on the preparation of 96 DNA samples.

**Normalization-free.** There is no need to control sample input or PCR cycles. The workflow auto-normalizes and produces similar amounts of sequencing reads across all samples regardless of different DNA inputs (Figure 2). Just pool by equal volumes!



**Figure 2.** The *Quick-16S™ Plus* NGS Library Prep Kit (V4) results in similar amounts of reads across different input amounts without normalization. 0.01-100 ng of ZymoBIOMICS™ Microbial Community DNA Standard and fecal DNA were used as inputs. Libraries were pooled by equal volumes (2 µl each) without further normalization and sequenced using the MiSeq® Reagent Kit v2 (300-cycle). The CV (coefficient of variation) is the ratio of the standard deviation to the mean with lower values corresponding to less dispersion around the mean.

**High Quality NGS Library.** The workflow has been optimized to minimize amplification bias (Figure 3).



**Figure 3.** Benchmarked performance with ZymoBIOMICS™ Microbial Community DNA Standard. Bacterial composition profiles are accurate with inputs from 100 ng down to 0.01 ng.

# Protocol

## Before Starting

- ✓ **Sample Quantity Requirement.** To ensure color balance in index sequencing, a minimum of 9 samples per run is recommended.
- ✓ **Input DNA Guidelines.** All DNA samples should be free of PCR inhibitors.<sup>1</sup> The 1-Step PCR reaction can accommodate DNA inputs of up to 100 ng but reducing inputs to  $\leq 10$  ng is recommended for optimal performance and robustness against potential PCR inhibition.
- ✓ **Sequencing Primer Requirement.** This kit utilizes custom sequencing primers that must be spiked into your Illumina cartridge prior to sequencing.

## Section 1: 1-Step PCR

1. Add 2  $\mu$ l of your DNA samples into the individual wells of the V4 Targeted Amplification Plate. Include a positive and negative control in the plate.
2. Apply an adhesive PCR plate seal. Mix the plate on a plate shaker and centrifuge in a plate spinner<sup>2</sup>.
3. Place plate in a real-time thermocycler<sup>3</sup> and run the program shown below:

Temperature	Time	
95°C	10 min	} 42 cycles <sup>4</sup>
95°C	30 sec	
55°C	30 sec	
72°C	3 min	
Plate read	-	
4°C	Hold	

*(Continued on next page.)*


<sup>1</sup> DNA that contains potent PCR inhibitors such as polyphenolics, humic/fulvic acids, tannins, melanin, etc. can be quickly cleaned using the [OneStep™ PCR Inhibitor Removal Kit](#).

<sup>2</sup> PCR reactions can be pipette mixed if a plate shaker is not available.

<sup>3</sup> A real-time thermocycler is recommended as it enables QC of the library prep of all wells. A non-quantitative system can be used if absolute quantification is not needed.

<sup>4</sup> The number of cycles can be adjusted to further reduce library prep time. See **Appendix B** for more details.

4. Monitor and QC the library preparation when running the reaction on a real-time thermocycler.<sup>1</sup>
  - a. For example, a sample that is expected to amplify and shows little or no amplification may indicate an error in the reaction setup (See the Troubleshooting Guide).
  - b. The negative control should not amplify before 35 cycles.<sup>2</sup> Earlier amplification of negative control may indicate process contaminations.
  - c. An example of qPCR amplification with controls is shown in Figure 5 below.



**Figure 5. qPCR Amplification Example with Positive and Negative Controls.** Serial dilutions of fecal DNA (black) from 100 ng to 100 fg were amplified on a Bio-Rad CFX96™ Real-Time PCR Detection System. The positive (blue) amplified at 16.59 and negative (red) amplified at 40.12. Baseline threshold was set at 1200 RFU.

5. Once the samples have cooled to 4°C, stop the program. Centrifuge plate in a plate spinner to collect condensation in wells and place plate on ice. Proceed to [Section 2](#), or store plate at ≤-20°C for later use. The plate will be stable at -20°C for up to a month.

## Section 2: Pooling by Equal Volume

Add 50 µl of **PCR Inactivation Solution** into a new microcentrifuge tube.<sup>3</sup> Pool equal volumes (2-8 µl) of PCR products from each well of the plate from [Section 1](#) into the tube<sup>4</sup> and mix well. If the total pooled library exceeds 500 µl, proceed to [Section 3](#) with a maximum of 500 µl of the well-mixed library.

<sup>1</sup> If real-time PCR was not used, after amplification perform PCR cleanup for a few samples plus positive control. Analyze on a TapeStation® to confirm correct amplicon size (~388 bp).

<sup>2</sup> The PCR program runs for 42 cycles, so it is normal to see some amplification from the negative control. The negative control should be sequenced together with other samples. If appropriate for your project, the taxa from the negative control can be subtracted from the analysis.

<sup>3</sup> If pooling from 384-well plate, a multi-channel pipette may be used to pool into 8-strip tubes pre-filled with 7 µl of **PCR Inactivation Solution** before consolidating into a microcentrifuge tube.

<sup>4</sup> There is no need for additional normalization procedures.



### Section 3: Final Library Clean-up

1. Equilibrate the **Select-a-Size™ MagBeads** to room temperature (15-30°C). Resuspend 30 µl of the **Select-a-Size™ Magbead Concentrate** in 1 ml of the **Select-a-Size™ Magbead Buffer**.<sup>1</sup> Resuspend the magnetic particles by vigorously shaking until homogenous.
2. Add **Select-a-Size™ MagBeads** to the pooled library from [Section 2](#) at a ratio of 0.8x volume. For example, add 400 µl of Select-a-Size™ MagBeads to 500 µl of the pooled library and PCR Inactivation Solution mixture.
3. Mix thoroughly by pipetting or vortexing until homogenous. Incubate for 5 minutes at room temperature.
4. Place the sample on a magnetic rack<sup>2</sup> and incubate for 3-10 minutes at room temperature, or until the magnetic beads have fully separated from solution.
5. Once the beads have cleared from solution, remove and discard the supernatant.<sup>3</sup>
6. While the beads are still on the magnetic rack, add 1 ml of **DNA Wash Buffer**. Remove and discard the supernatant. Repeat this step.
7. While the beads are still on the magnetic rack, aspirate out any residual buffer with a 10 µl pipette tip.
8. Remove tube from the magnetic rack and keep the cap open for 3 minutes at room temperature to dry the beads.
9. Add 10-100 µl<sup>4</sup> of **ZymoBIOMICS™ DNase/RNase Free Water** to the beads and pipette mix thoroughly. Incubate at room temperature for 2 minutes.
10. Place the sample on a magnetic rack and incubate for 1 minute at room temperature, or until the magnetic beads have fully separated from eluate.
11. Transfer supernatant to a clean microcentrifuge tube. Proceed to [Section 4](#).

<sup>1</sup> Once the concentrate and buffer have been mixed, the mix can be stored at 4°C for a maximum of 3 months.

<sup>2</sup> Alternatively, the provided Magnetic Rods can be used.

<sup>3</sup> Avoid aspirating any beads when removing the supernatant. To best prevent this, leave 2-5 µl of liquid behind.

<sup>4</sup> If pooling fewer than 10 samples, use 10 µl for elution.

## Section 4: Library Quantification

Use a fluorescence-based method (Qubit® dsDNA HS Assay Kit recommended) to quantify the final library. Using a final amplicon size of 388 bp, convert ng/μl to nM using the equation below.

$$\frac{\text{concentration in ng/}\mu\text{l}}{660 \text{ g/mol} \times \text{average library size in bp}} \times 10^6 = \text{concentration in nM}$$

If preferred, a qPCR-based method for quantification such as the KAPA® Library Quantification kit may be used.

### **DNA Fragment Analysis (Not Required)**

If a fragment analyzer (e.g., TapeStation®) is used to analyze the final library, there may be a lack of a tight band at ~388 bp. Because of the library prep design, some library products have run through additional PCR cycles and might not anneal well. This is perfectly fine for sequencing. Double-stranded DNA will be denatured into single strands before loading onto the sequencer.

### **This is your final 16S library**

The ultra-pure pooled library DNA is now ready for use or storage at ≤-20°C. Refer to platform-specific guidelines for preparation for sequencing.

### **Illumina MiSeq® Setup:**

The MiSeq® Reagent Kit v2 (300 Cycle) with 15% PhiX spike-in is recommended. See Appendices D for assistance with sample sheet setup. Remember to set the index size to 10 bp. Proceed to [Section](#)

[5.](#)

## Section 5: Denaturation and Dilution

Using the recommended loading concentrations below, denature and dilute your final library according to the Illumina® Denature and Dilute Protocol for your specific sequencing platform. A 15% PhiX spike-in is recommended across all Illumina® Platforms.

Recommended final library loading concentrations:

Illumina® Platform	Loading Concentration
MiSeq v2	4 pM
MiSeq v3	8 pM
MiniSeq	1 pM
NovaSeq	150 pM

The final loading concentration may need to be further optimized for your specific machine to achieve optimal cluster densities.

Proceed to [Section 6](#).

## Section 6: Spiking Custom Primers into Illumina Reagent Cartridge

1. Refer to **Appendix E** for the primer positions.
2. Using a pipette tip, pierce the foil at the appropriate position on the cartridge.
3. Using a Pasteur pipette, aspirate the contents of the cartridge well and transfer into an empty microcentrifuge tube.
4. Add the appropriate volume of 100  $\mu$ M custom primer into the tube. Vortex and quick spin.
5. Using a P1000, transfer the contents of the tube back to its original position on the cartridge.
6. Repeat for all the necessary primers.
7. Load the final library into the proper position.
8. Look at the bottom of the cartridge and make sure that there are no bubbles in the reagent wells.
9. Your cartridge is now ready for sequencing.

## **Appendix A: Removal of PCR Inhibitors from Starting DNA**

The input DNA samples for the **Quick-16S™ Plus NGS Library Prep Kit (V4)** must be free of PCR inhibitors such as polyphenolics, humic/fulvic acids, tannins, melanin, etc. To further remove PCR inhibitors from purified DNA samples, Zymo Research recommends performing a one-step cleanup with the **OneStep™ PCR Inhibitor Removal Kit**. Additional information can be found by visiting the **D6030** Product Page at [www.zymoresearch.com](http://www.zymoresearch.com).

## **Appendix B: Adjusting PCR Cycles**

If samples, excluding the no template control, consistently plateau before 42 cycles, the PCR cycles can be adjusted to further reduce library prep time. Using previous data that is representative of your samples, determine the Ct which all samples, excluding the no template control, plateau. If samples plateau before 25 cycles, set the PCR program to 25 cycles. If the samples plateau before 30 cycles, set the PCR program to 30 cycles. If the samples plateau before 35 cycles, set the PCR program to 35 cycles.

## **Appendix C: Sequencing Recommendations for Other Platforms**

For sequencing platforms not listed, we recommend diluting your final library to a concentration lower than the recommended concentration in the standard Illumina protocol for your specific platform to prevent over clustering. Increase the loading concentration until you achieve optimal cluster densities. We recommend 15% PhiX spike-in across all platforms and reagent kits. Final library concentrations and PhiX spike-in percentages may be further adjusted to uniquely optimize cluster density and sequencing read quality respectively.

## **Appendix D: Illumina MiSeq® Sample Sheet Setup**

A template for the Illumina MiSeq® sample sheet is available for download by visiting the Documentation section of the Product Page at [www.zymoresearch.com](http://www.zymoresearch.com). Fill in the project and sample information in the highlighted fields, then save the file in comma-separated values (CSV) format for use with the Illumina MiSeq®.

## Appendix E: Illumina Sequencing Primer Positions


Sequencing Kit	Illumina Primer (name)	Cartridge Position	Vol of Custom Primer Spike In ( $\mu$ l)
MiniSeq High Output (300 cycles)	Read 1	24	1.65
	Read 2	25	1.83
	Index 1 and Index 2	28	2.46
MiniSeq Mid Output (300 cycles)	Read 1	24	1.65
	Read 2	25	1.83
	Index 1 and Index 2	28	2.46
MiSeq v2 and v3	Read 1	12	3.4
	Index 1	13	3.4
	Read 2	14	3.4
NextSeq 500/550 High Output (300 Cycles)	Read 1	20	5.19
	Read 2	21	5.94
	Index 1 and Index 2	22	8.49
NextSeq 500/550 Mid Output (300 Cycles)	Read 1	20	3.99
	Read 2	21	4.56
	Index 1 and Index 2	22	6.27
NovaSeq v1.0 SP (300 Cycles)	Read 1	24	12
	Index 1	23	15
	Read 2	13	6
NovaSeq v1.0 S1 and S2 (300 Cycles)	Read 1	24	12
	Index 1	23	15
	Read 2	13	6
NovaSeq v1.0 S4 (300 Cycles)	Read 1	24	21.9
	Index 1	23	15
	Read 2	13	10.5
NovaSeq v1.5 SP (300 Cycles)	Read 1	24	12
	Index 1 and Index 2	23	15
	Read 2	13	6
NovaSeq v1.5 S1 and S2 (300 Cycles)	Read 1	24	12
	Index 1 and Index 2	23	15
	Read 2	13	6
NovaSeq v1.5 S4 (300 Cycles)	Read 1	24	21.9
	Index 1 and Index 2	23	15
	Read 2	13	10.5

\*For reagent kits that list both Index 1 and Index 2, the primers should be spiked into the same position.

For example: For the MiniSeq High Output Kit, 2.46  $\mu$ l of Index 1 Sequencing Primer and 2.46  $\mu$ l of Index 2 Sequencing Primer should be spiked into cartridge position 28.

# Explore Other Microbiome Products

- ✓ To collect and transport samples at ambient temperatures:




DNA/RNA Shield™ and Collection Devices	
1X Reagent #R1100	For sample lysis and stabilization of DNA/RNA
2X Concentrate #R1200	Reagent concentrate (2X) for use with liquids at 1:1 ratio
Fecal Collection Tube #R1101	15 mL container (prefilled with 9 mL DNA/RNA Shield™). Direct collection of up to 1g or 1 mL stool
Collection Tube w/ Swab #R1106	12 x 80 mm screwcap container filled with 1 mL DNA/RNA Shield™ and sterile swab for specimen collection

- ✓ Unbiased and inhibitor-free DNA and RNA extraction (high-throughput and automatable) for microbial profiling:



ZymoBIOMICS™ DNA and RNA Kits	
DNA Miniprep #D4300	Up to 25 µg DNA
DNA Microprep #D4301	Up to 5 µg DNA
MagBead DNA #D4302	Automatable (Tecan, Hamilton, Kingfisher, etc.)
96-Well DNA #D4309	Spin-plate
DNA/RNA Miniprep Kit #R2002	Up to 100 µg DNA/RNA

- ✓ Microbial standards and references for profiling quality control, benchmarking, positive controls, and to assess performance of entire microbiomic/metagenomic workflows:



ZymoBIOMICS™ Standards and Reference Materials	
Microbial Community Standard #D6300	Contains 8 bacteria and 2 yeasts for QC and method optimization
Microbial Community DNA Standard #D6305	Contains 8 bacteria and 2 yeasts DNA for bioinformatics optimization
Gut Microbiome Standard #D6331	Contains 21 different human gut strains for method benchmarking
Fecal Reference with TruMatrix™ Technology #D6323	Contains real human fecal material for benchmarking and improved data reproducibility

# Troubleshooting Guide

Problem	Possible Causes and Suggested Solutions
<b>Background Contamination</b>	<p>Workspace contamination:</p> <ul style="list-style-type: none"> <li>- Clean workspace, microcentrifuge, and pipettes with 10% bleach routinely to avoid contamination.</li> <li>- Use of kit in exposed environment without proper filtration can lead to background contamination. Check pipettes, pipette tips, microcentrifuge tubes, workspace, etc. for contamination</li> <li>- Make sure all reagent tubes and bottles are properly sealed for storage. Use of these outside a clean room or hood can result in contamination.</li> </ul>
<b>Loss of Volume during PCR</b>	<p>Adhesive seal:</p> <ul style="list-style-type: none"> <li>- A loosened adhesive seal on the PCR plate can lead to sample evaporation. Ensure that the plate seal is secure on every well during targeted sequence amplification.</li> </ul> <p>Lid pressure:</p> <ul style="list-style-type: none"> <li>- Inconsistent lid pressure. Ensure that the lid pressure on the real-time quantitative PCR instrument is consistent over the PCR plate according to the manufacturer's recommendation.</li> </ul>
<b>Unexpected or No Amplification of DNA Sample During PCR Program in <a href="#">Section 1</a></b>	<p>Sample with high microbial DNA concentration:</p> <ul style="list-style-type: none"> <li>- Reaction setup error. A sample that is expected to amplify but shows little or no amplification during the PCR program in <a href="#">Section 1</a> may indicate an error in the reaction setup. Use a new aliquot of the sample and repeat <a href="#">Section 1</a>.</li> <li>- Sample may contain high levels of PCR inhibitors. See Appendix A on how to remove these and repeat <a href="#">Section 1</a>. Additionally, samples may be diluted to lower concentration (&lt;10 ng/μl) to see if that improves amplification efficiency.</li> </ul> <p>Sample with low microbial DNA concentration:</p> <ul style="list-style-type: none"> <li>- Check negative control. A sample with little microbial DNA may not amplify before the negative control. Either use more concentrated DNA or use more DNA volume during reaction setup.</li> </ul> <p>Abnormal qPCR curves:</p> <ul style="list-style-type: none"> <li>- Proceed as normal. Abnormal qPCR amplification curves may occur, and this is normal performance. This is usually a slight dip in RFU (forming a small "hump").</li> </ul>
<b>Diminished Amplicon Bands in Library Analysis</b>	<p>No single amplicon peak and/or high background:</p> <ul style="list-style-type: none"> <li>- Proceed as normal. There may be a lack of a single band and/or high background if using TapeStation® or similar methods to determine amplicon size. This is normal and is part of the library prep design. Do not use this sizing and quantification data. To properly quantify library, use a fluorescence-based method and calculation in <a href="#">Section 4</a>.</li> </ul>

For technical assistance, please contact 1-888-882-9682 or email [tech@zymoresearch.com](mailto:tech@zymoresearch.com)







# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---



## **100% satisfaction guarantee on all Zymo Research products, or your money back.**

Zymo Research is committed to simplifying your research with quality products and services. If you are dissatisfied with this product for any reason, please call 1(888) 882-9682.

---

Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

<sup>™</sup> Trademarks of Zymo Research Corporation, <sup>®</sup> Registered Trademarks of Zymo Research Corporation  
Other trademarks: Bio-Rad CFX96<sup>™</sup> Real-Time PCR Detection System (Bio-Rad Laboratories, Inc.), MiSeq<sup>®</sup> (Illumina, Inc.), TapeStation<sup>®</sup> (Agilent Technologies, Inc.), Qubit<sup>®</sup> (Thermo Fisher Scientific), Kapa<sup>®</sup> (Roche), TruMatrix<sup>™</sup> (The BioCollective).



ZYMO RESEARCH

*The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®*