



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

Quick-DNA/RNA™ Pathogen MagBead

Catalog Nos. R2145 & R2146

Highlights

- High-throughput, magnetic-bead based purification of pathogen (virus, bacteria, protozoa) DNA/RNA from a wide variety of vectors (mosquitoes, fleas, ticks, *etc.*) and tissue types (mammals, birds, *etc.*) stored in DNA/RNA Shield™
- DNA/RNA is ready for any sensitive downstream applications (e.g., Next-Gen sequencing, RT/PCR, *etc.*)

Contents

Product Contents	1
Specifications	1
Product Description	2
Buffer Preparation	3
Protocol	3-4
Appendix	5
Ordering Information	6

Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product please contact us.

Note: Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Product Contents

Quick-DNA/RNA™ Pathogen MagBead (Kit Size)	R2145 (96 preps)	R2146 (4x 96 preps)
DNA/RNA Shield™	2x 50 ml	2x 250 ml
Pathogen DNA/RNA Buffer¹	100 ml	4x 100 ml
Proteinase K w/ Storage Buffer²	5 mg	20 mg
MagBinding Beads	3 ml	12 ml
MagBead DNA/RNA Wash 1 (concentrate)³	30 ml	120 ml
MagBead DNA/RNA Wash 2 (concentrate)⁴	20 ml	80 ml
DNase/RNase-Free Water	10 ml	2x 30 ml
Instruction Manual	1	1
ZR Bashing Bead™ Lysis Tubes (sold separately)	S6014-50 (1x 50 pack or 8x 50 pack)	

Storage Temperature – Store all kit components (*i.e.*, buffers, columns) at room temperature.

¹ Add beta-mercaptoethanol to 0.5% (v/v) *i.e.*, add 500 µl or 1 ml per 100 ml or 200 ml **Pathogen DNA/RNA Buffer**, respectively.

² Add 260 µl or 1,040 µl **Proteinase K Storage Buffer** per vial to reconstitute the lyophilized **Proteinase K**, 5 mg or 20 mg respectively.

Vortex to dissolve and store frozen aliquots.

³ Add 20 mL (R2145) or 80 ml (R2146) of isopropanol to the **MagBead DNA/RNA Wash 1** concentrate.

⁴ Add 30 mL (R2145) or 120 ml (R2146) of isopropanol to the **MagBead DNA/RNA Wash 2** concentrate.

Specifications

- **Sample Sources** – vectors (mosquitoes, fleas, ticks, other tough-to-lyse insects) and tissue types (animal tissue, plants, other hosts) processed and stored in DNA/RNA Shield™.
- **Sample Size** – ≤200 µl (or ≤10 mg)
- **Format** – 10 µg DNA/RNA binding capacity (per 20 µl MagBinding Beads), ≥50 µl elution volume
- **Size Range** – 50 nt to ~200 kb
- **Purity** – DNA/RNA is ready for any sensitive downstream applications (e.g., Next-Gen sequencing, RT/PCR, *etc.*)
- **Materials** (sold separately) –
 - ZR BashingBead Lysis Tubes (S6014; 0.1/2.0 mm)
 - ZR-96 MagStand (P1005)
 - Collection Plate (C2002; capacity 1.2 ml/well)
 - 96-Well Block (P1001; capacity 2 ml/well)
 - Elution Plate (C2003; capacity 0.35 ml/well)
 - Cover Foil (C2007)

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com


Product Description

Quick-DNA/RNA™ Pathogen MagBead kit is designed for high-throughput purification of pathogen (virus, bacteria, protozoa) DNA and RNA from a wide variety of vectors (mosquitoes, fleas, ticks, etc.) and tissue types (mammals, birds, etc.) collected, transported and stored in **DNA/RNA Shield™**.


The kit features a storage/lysis buffer system and can be combined with high density ZR BashingBead™ Lysis Tubes (optional) to facilitate complete homogenization of hard-to-lyse samples for efficient nucleic acid isolation. Small (>50 nt) and large (>200 kb) DNA and RNA are bound to magnetic beads, washed and then eluted.

The isolated high-quality nucleic acids are suitable for all downstream applications such as Next-Gen sequencing, hybridization-based and RT/PCR detection.


For Assistance, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.



Inhibitor-free detection of West Nile Virus in mosquitoes



High-sensitivity detection of HIV-1 virus at low titer



(A) West Nile Virus (spiked-in mosquito homogenate), (B) HIV-1 viral RNA particles (spiked-in plasma), purified using the **Quick-DNA/RNA™ Pathogen** kit and detected by RT-qPCR.

*optional
ZR BashingBead Lysis Tubes 2.0 mm + 0.1 mm (catalog no. S6014)



All reagents and steps should be performed at room temperature, unless specified.

Buffer Preparation

- Add 500 μ l or 1 ml beta-mercaptoethanol (user supplied) per 100 ml or 200 ml **Pathogen DNA/RNA Buffer**, respectively (final concentration 0.5% (v/v)).
- Add 20 mL (R2145) or 80 ml (R2146) of isopropanol to the **MagBead DNA/RNA Wash 1** concentrate.
- Add 30 mL (R2145) or 120 ml (R2146) of isopropanol to the **MagBead DNA/RNA Wash 2** concentrate.
- Add 260 μ l or 1,040 μ l **Proteinase K Storage Buffer** per vial to reconstitute the lyophilized **Proteinase K**, 5 mg or 20 mg. Vortex to dissolve and store frozen aliquots.

Protocol

- All steps should be performed at room temperature, unless specified.

Sample Preparation

1. Add **DNA/RNA Shield™** to the sample as recommended (table below) in a nuclease-free tube (not provided) and mix well:

		add DNA/RNA Shield™ #R1100
Insects ¹ (tough-to-lyse, ticks, mosquitoes, fleas, deer flies, etc.)	≤ 10 mg	800 μ l
Tissue (mammals, birds, plants)		

For biological samples² collected and stored in DNA/RNA Shield™ (i.e., #R1200, R1101, R1150, R1104, R1107, R1109, etc.), proceed to step 3 below.

2. Optional: Mechanical homogenization with the ZR BashingBead™ Lysis Tube³ and a high-speed cell disruptor⁴ is recommended for tough-to-lyse insects, animal tissue, plants, etc.

	homogenization time (high-speed)
insects (tough-to-lyse, ticks, mosquitoes, fleas, deer flies, etc.)	3-5 minutes
tissue (mammals, birds, plants)	30-60 seconds

3. To remove particulate debris or precipitation, centrifuge at 10,000-16,000 x g for 1 minute. Transfer up to 200 μ l of the cleared supernatant⁵ into a nuclease-free tube or well/plate (not provided).

Notes:

¹ See Appendices, page 5 for specific number of insect input.

² Plasma, serum, whole-blood, urine, fecal, swab, saliva, cell suspension, culture media, etc.

³ S6014-50, S6014-200

⁴ Required homogenization time will vary depending on the device and application. For high-speed cell disruptors (e.g., FastPrep® - 24, TerraLyzer™ Sample Processor or similar), samples can be processed in ≤ 5 minutes. For low-speed cell disruptors (e.g., Disruptor Genie™, or standard benchtop vortexes), processing can be ≤ 20 minutes long. See manufacturer's literature for operating information.

⁵ Up to 200 μ l liquid sample can be processed per prep.

DNA/RNA Purification

4. Add 2 μ l **Proteinase K** and mix well.
 5. Add 400 μ l **Pathogen DNA/RNA Buffer**⁶ to each 200 μ l sample⁷ and mix well⁸.
 6. Add 20 μ l **MagBinding Beads** and mix well⁸ for 10 minutes. Important: MagBinding Beads settle quickly, ensure that beads are kept in suspension while dispensing.
 7. Transfer the plate to a magnetic stand until beads have pelleted, then aspirate⁹ and discard the cleared supernatant.
 8. Add 500 μ l **MagBead DNA/RNA Wash Buffer 1** and mix well⁸.
 9. Transfer the plate to a magnetic stand until beads have pelleted, then aspirate⁹ and discard the cleared supernatant.
 10. Add 500 μ l **MagBead DNA/RNA Wash Buffer 2** and mix well⁸.
 11. Transfer the plate to a magnetic stand until beads have pelleted, then aspirate⁹ and discard the cleared supernatant.
 12. Add 500 μ l ethanol (95-100%) and mix well⁸.
 13. Transfer the plate to a magnetic stand until beads have pelleted, then aspirate⁹ and discard the cleared supernatant.
 14. Repeat steps 12 and 13.
 15. Dry the beads at room temperature for 10 minutes or until fully dry¹⁰.
 16. To elute DNA/RNA from the beads, add 50 μ l **DNase/RNase-Free Water** and mix well⁸.
- Alternatively, for highly concentrated DNA/RNA use ≥ 30 μ l volume.
17. Transfer the plate to a magnetic stand until beads have pelleted, then aspirate⁹ and dispense the eluted DNA/RNA to an elution plate.

The eluted DNA/RNA¹¹ can be used immediately or stored frozen.

All reagents and steps should be performed at room temperature, unless specified.

Notes:

⁶ To ensure efficient lysis and deproteinization, up to 5 volumes of Pathogen DNA/RNA Buffer can be used per 200 μ l liquid sample.

⁷ Up to 200 μ l liquid sample can be processed per prep.

⁸ For all buffer additions and to ensure beads are properly in suspension, **mix well** by pipetting up and down several times and/or by shaking (vortexing) at $\sim 1,300$ rpm.

⁹ Some beads will adhere to the sides of the well. When removing the supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered.

¹⁰ Beads will change in appearance from glossy black when still wet to a dull brown when fully dry. Alternatively, a heat block can be used (25-55°C).

¹¹ It is recommended to titrate the DNA/RNA eluate for downstream applications (i.e., RT/PCR, etc.).

Appendix

Insect Samples: Recommended Input

Sample type	Maximum specimen input per 800 μ l DNA/RNA Shield™
mosquito	≤ 50
tick	1 engorged of any species ≤ 5 flat adults ≤ 20 nymphs
flea	≤ 10
deer fly	1 adult

Automation Scripts

The **Quick-DNA/RNA™ Pathogen MagBead** kit is compatible with any automated platform. For automation scripts and related technical support, email automation@zymoresearch.com. In the subject line, please include “Automation Scripts”, instrument used and the product catalog number.

[Automation Script for MagMAX Express 24 – AM1836 protocol](#)

Ordering Information

Product Description	Catalog No.	Kit Size
Quick-DNA/RNA™ Pathogen MagBead	R2145	1x 96 Preps
	R2146	4x 96 Preps
Quick-DNA/RNA™ Pathogen MiniPrep	R1042	50 Preps
	R1043	200 Preps

For Individual Sale	Catalog No.	Amount
DNA/RNA Shield™	R1100-50	50 ml
	R1100-250	250 ml
Pathogen DNA/RNA Buffer	R1042-1-50	50 ml
	R1042-1-100	100 ml
Proteinase K w/ Storage Buffer Set	D3001-2-5	5 mg
	D3001-2-20	20 mg
MagBinding Beads	D4100-2-3	3 ml
	D4100-2-12	12 ml
MagBead DNA/RNA Wash 1 (concentrate)	R2130-1-30	30 ml
	R2130-1-120	120 ml
MagBead DNA/RNA Wash 2 (concentrate)	R2130-2-20	20 ml
	R2130-2-80	80 ml
DNase/RNase-Free Water	W1001-30	30 ml
	W1001-100	100 ml
	W1001-200	200 ml
ZR BashingBead™ Lysis Tubes (0.1 mm & 2.0 mm)	S6014-50	50 pack

ZYMO RESEARCH CORP.

DNA & RNA MADE SIMPLE



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com