

LiveReceptor™ GABA_AR <GABA_AR Labeling Reagent>

Catalog NO. FDV-0018B

Research use only, not for human or animal therapeutic or diagnostic use.

Product Background

Neurotransmitter receptors including glutamate receptors and GABA receptors etc. located on post-synapse in neuronal cells play various roles in brain functions. To understand physiological roles of neurotransmitter receptors, live cell imaging is one of the powerful approaches. Conventional imaging methods on live cells rely on a genetically engineered proteins fused with fluorescent proteins such as GFP. However, one serious problem is that the functions and movement of over-expressed neurotransmitter receptors with non-physiological tags are not precisely correlated with endogenous native receptors. The labeling methods for endogenous receptors are desirable to observe physiological functions of receptors.

LiveReceptor™ is the world first reagent series for target-specific receptor labeling. The principle of LiveReceptor™ is based on ligand-directed acylyl imidazole (LDAI) chemistry (ref.1,2). LDAI-based chemical labelling is driven by selective ligand-protein recognition, which facilitates an acyl substitution reaction of labeling reagents on nucleophilic amino acid residues including Lys, Ser and Tyr located near ligand-binding domain. After wash out, the labelled receptors which have free ligand-binding pockets are observed on live cells. Furthermore, based on pH-dependent fluorescent property of fluorescein, fluorescent signal of labeled receptors in endocytosis pathway are highly quenched and only cell surface receptors can be observed. LiveReceptor™ series are powerful tools to monitor reduction of cell surface receptors by endocytosis upon extracellular stimulation.




Figure 1. Principle of ligand-directed protein labeling

“ **LiveReceptor™ GABA_AR**” is a specific labeling reagent for cell-surface ion channel-type GABA receptor (GABA_AR) which is the key component for inhibitory synaptic regulation (ref.3). LiveReceptor™ GABA_AR has three domains including gabazine as an affinity ligand for GABA_AR, fluorescein and acyl imidazole. Only when gabazine binds to GABA_ARs, nucleophilic amino acid residues (Lys, Ser or Tyr) located near ligand-binding domain on GABA_ARs are attacked acyl imidazole and fluorescein is transferred into GABA_ARs. After removing excess reagents and resultant ligand moiety, labeled GABA_ARs can be observed in both live and fixed cells. The protocol is very simple, no genetic manipulation and additional treatment are required. Because LiveReceptor™ GABA_AR shows no cell membrane permeability, only cell surface GABA_ARs are labelled. Ref.3 indicates fluorescein-labeled GABA_AR by LiveReceptor™ has little effects on its ion channel capability.




Figure 2. Chemical structure of LiveReceptor™ GABA_AR

Description

Catalog Number: FDV-0018B

Size : 10 µg

Formulation : C₅₅H₅₁N₇O₁₃

Molecular weight : 1017.35 g/mol

Visibility : Orange lyophilized powder

Solubility : Soluble in DMSO

*This compound has water-solubility but it can be easily degraded in water and culture medium.

Please avoid store in the water.

Spectrum

Excitation/ Emission: 495/515 nm

*Compatible with FITC filter

Application

- Live cell imaging
- Immunocytochemistry with specific antibodies
- Immunoprecipitation with anti-fluorescein antibody
- Immunoblotting with anti-fluorescein antibody
- Drug screening for competitive GABA_AR antagonist (GABA-site)

Reconstitution and Storage

Reconstitution :

Reconstitute at 0.1 mM (x100) - 1 mM (x1000) in 100% DMSO. Please optimize the final concentration of DMSO depended on your experiments. Before reconstitution, please spin down to collect the orange lyophilized powder on the bottom of a tube. Carefully add DMSO into the tube and vigorously mix to completely dissolve the powder.

Storage:

(powder) Store at -20°C. Protected from light.

(solution) DMSO stock solution is stable at least for 1 year at -80°C. Please make aliquots and avoid freeze and thaw. Protected from light.

How to use


General procedure for GABA_AR labeling

1. Prepare 1 μM of LiveReceptor™ GABA_AR in the appropriate medium
* Note : Serum-free media are highly recommended. This compound is not stable in the medium. Please prepare assay solution at time of use.
2. Replace media of cultured cells to LiveReceptor™ containing medium.
3. Culture cells with LiveReceptor™ GABA_AR for 1-4 hours at 17-37°C
4. After labeling, wash cells several times or perfused continuously to remove excess reagents.
5. Labelled GABA_ARs can be observed.

Application data


Live cell imaging of labelled GABA_ARs in GABA_AR-expressed HEK293

GABA_AR (α1/β3/γ2)-expressed HEK293 cells were treated with 1 μM of LiveReceptor™ GABA_AR in the absence or presence of 100 μM gabazine, a GABA_AR selective inhibitor, for 3 hour at 37°C and washed out three times with the basal medium. (scale bars, 20 μm)



Validation of GABA_AR labeling by western blotting

GABA_AR (α1/β3/γ2)-expressed HEK293 cells were treated with 0.01-1 μM of LiveReceptor™ GABA_AR in the absence or presence of 100 μM gabazine, a GABA_AR selective inhibitor, for 3 hour at 37°C. After wash cells, the cells were lysed and proteins were applied into SDS-PAGE and western blotting using anti-fluorescein (FL) antibody or anti-GABA_AR α1 subunit.



Reference

1. Fujishima *et al.*, *J. Am. Chem. Soc.*, **134**, 3961-3964 (2012). Ligand-directed acyl imidazole chemistry for labeling of membrane-bound proteins on live cells.
2. Miki *et al.*, *Chem. Biol.*, **21**, 1013-1022 (2014). LDAI-based Chemical Labeling of Intact Membrane Proteins and its Pulse-Chase Analysis under Live Cell Conditions.
3. Yamaura *et al.*, *Nat. Chem. Biol.*, **12**, 822-830 (2016). Discovery of allosteric modulators for GABA_A receptors by ligand-directed chemistry.

Disclaimer/免責事項

This product has been commercialized by Funakoshi Co., Ltd. based on the results of academic research, and the advertisement text, figures and manuals (hereinafter “Product information”) have been prepared based on published research reports on January, 2019. The academic interpretation at the time of creation of the Product Information may change in accordance with future developments in the relevant research field and expansion of various scientific findings, and the latest version and certainty of the Product Information are not guaranteed. The specifications of this product and the Product Information are subject to change without notice. Please contact us for the latest information.

本商品は学術研究成果を基にフナコシ株式会社が商品化したもので、2019年1月時点における公開研究報告を基に広告文章およびマニュアル(以下、商品資料)を作成しています。今後の当該研究分野の発展および各種学術知見の拡大にともない、商品資料作成時の学術的解釈が変更になる可能性があり、最新性・確実性を保証するものではありません。また、本商品の仕様および商品資料を予告なく変更する場合がございます。最新の情報に関しましては、弊社までご確認いただけますようお願い申し上げます。



E-mail Newsletter
Sign Up

Japanese



English



Related product

LiveReceptor™ AMPAR <Endogenous AMPAR Labeling Reagent>

LiveReceptor™ AMPAR is a specific labeling reagent for AMPA-type glutamate receptor, AMPAR. Live imaging of cultured neuron and slice tissues were validated.


Catalog No. FDV-0018A

Size 10 µg

Data examples

- Live cell imaging of labelled endogenous AMPARs in cultured neurons

Cultured hippocampal neurons were treated with 1 µM of LiveReceptor™ AMPAR (in left) or Fluorescein-conjugated PFQX as negative control (in right) for 1 hour at 17°C and washed out three times with the basal medium. Dendritic spin-like punctual structures were observed on live cells by specifically LiveReceptor™. (scal bars, 10 µm)



LiveReceptor™ mGluR1 <Endogenous mGluR1 Labeling Reagent>

LiveReceptor™ mGluR1 is a specific labeling reagent for metabotropic glutamate receptor 1, mGluR1. Live imaging of slice cerebellum tissue was validated.


Catalog No. FDV-0018C

Size 10 µg

Data examples

- Live cell imaging of labelled mGluR1 in acute mouse cerebellum brain slice

Acute mouse cerebellum brain slices from 3 week-old mouse treated with 10 nM LiveReceptor™ mGluR1 in artificial cerebrospinal fluid (ACSF) for 4 hours at RT. After then, slices were washed three times by ACSF and fluorescent signal was observed by epi-fluorescent microscopy. Strong fluorescent signal was observed in molecular layer and Purkinje cells. When the slice was pretreated with FITM, an allosteric mGluR1 inhibitor, fluorescent signal was clearly suppressed.



Featured products

LipiORDER™ <Membrane Lipid Order Imaging Dye>


LipiORDER™ is a solvatochromic dye for membrane lipid order imaging. LipiORDER™ exhibits green fluorescence with Lo phase and exhibits red fluorescence with Ld phase. The ratiometric analysis (F_{red}/F_{green}) enables the quantitative visualization of membrane lipid order.

Catalog No. FDV-0041

Size 0.1 mg

Features

- Recommended Ex/Em: ~405 nm / 500-550 nm (Green channel) and 550-650 nm (Red channel)
- Enable to quantitatively monitor lipid order on plasma and inner membranes in live cells
- Highly photostable and cellularly stable compared with similar conventional dyes.



LipiDye™ II <Lipid Droplet Live Imaging>


LipiDye™ II is a highly sensitive lipid droplet staining dye with extremely photostable property. This dye is the second generation of our previous reagent, LipiDye™. This dye allows us to detect small lipid droplets (<1 μm) in non-adipocytes and to apply into long-term live cell imaging for dynamic lipid droplet movements.

Catalog No. FDV-0027

Size 0.1 mg

Features

- Recommended Ex/Em: 400-500 nm / 490-550 nm
- Enable to detect <1 μm lipid droplets
- Suitable for long-term live cell imaging
- Extremely photostable compared with conventional dyes
- Compatible with both live and fixed cells



FAOBlue™ <Fatty Acid Oxidation Detection Reagent>


FAOBlue™ is a cell-based fatty acid beta-oxidation (FAO) detection dye which emits blue fluorescence upon cellular FAO activity.

Catalog No. FDV-0033

Size 0.2 mg

Features

- Ex/Em: ~405 nm / 460 nm
- Enable to directly detect cellular FAO activity in live cells
- Apply quantitative comparison of FAO activity between different cell types
- Can monitor the drug-induced change of FAO activity



 **funakoshi**
FRONTIERS IN LIFE SCIENCE

URL: <http://funakoshi.co.jp>
9-7 Hongo 2-Chome, Bunkyo-ku, Tokyo 113-0033