

PolyamineRED <Intracellular Polyamine Detection Reagent>

Catalog NO. FDV-0020

Research use only, not for human or animal therapeutic or diagnostic use.

This product has been commercialized with the support of Biofunctional Synthetic Chemistry Laboratory, Cluster for Pioneering Research, RIKEN

Product Background

The polyamine species (Figure 1), including putrescine, spermidine and spermine etc. and its acetyl derivatives, are one of the essential class of metabolites which have linear alkyl structure with two or more amines. Polyamines are found in all living organisms with high concentration, from sub-millimolar to millimolar, in the cells. Polyamines have polycationic properties and shows an enormous number of biological functions. For example, polyamines interact with DNA/RNA in the nuclear and regulate gene expression. Polyamines also interact with negatively charged proteins and control its function. The major source of polyamines is an amino acid ornithine. In the case of mammalian, ornithine is converted to putrescine by ornithine decarboxylase (ODC), followed by synthesizing spermidine and spermine. Because ODC is highly expressed in cancer cells, polyamines are considered as one of the cancer marker. Several detection methods of polyamines are developed to date but most of the methods are commonly low-throughput systems using HPLC with polyamine standard compounds. To clear biological functions of polyamines in the cells, the cell-based assay with easy- and high-throughput-procedures is desired.

PolyamineRED is the world first reagent for detecting intracellular polyamines without any pre-treatment and cell lysis. PolyamineRED is a TAMRA (tetramethylrhodamine)-conjugated derivative of glycine propargyl ester which specifically reacts with linear primary alkylamine but not react with secondary amines, bulky amines including amino acids nor monoamines. PolyamineRED has cell-penetrating properties, specifically reacts with polyamines inside the cells and labelled polyamines with red fluorescent dye TAMRA. (Figure 2).




Figure 1. Major polyamine species

Description

Catalog Number: FDV-0020
Size : 0.5 mg
Molecular weight : 611 g/mol
Solubility : Soluble in DMSO
Fluorophore : TAMRA
(red fluorescent dye)
Ex/Em: 560 nm/585 nm
*TAMRA filter sets are available.




Figure 2. Principle of PolyamineRED

Reconstitution and Storage

Reconstitution : stock solution in 100% DMSO.

Storage (solution) :

Store powder at -20°C.

After reconstitution in DMSO, aliquot and store at -20 °C. Avoid repeated freeze-thaw cycles.


Protect from light.

How to use

General procedure of detection of intracellular polyamines

1. Prepare 10-30 μM PolyamineRED in fresh medium
 2. Remove culture medium, wash cells by PBS twice and add PolyamineRED-containing medium to cells
 3. Culture cells for at least 10 min
 4. Wash cells with PBS 3 times
 5. Fixed cells with paraformaldehyde
- Note: MeOH-fixation is not available. Please fix cells by formaldehyde.
6. Additional staining such as DAPI staining or immunocytochemistry with antibodies of interest are available.

1) Preparation of PolyamineRED-containing medium



Reference data

Selectivity of glycine propagyl ester to polyamines

Benzoyloxycarbonyl glycine propagyl ester as a model molecule was selectively reacted with polyamines. Reactant of epinephrine, an example of monoamine, and lysine, an example of amino acid, were rarely detected. Reactivity for polyamines depends on the length of polyamine and double linkage products were observed from spermine (4 amino groups) and spermidine (3 amino groups).

Table Selectivity of glycine propagyl ester to biological amines


amine	Reaction products			Hydrolysis product	Non reacted product
	Total	Single linkage	Double Linkage		
Spermine	82%	59%	23%	17%	1%
Spermidine	78%	67%	11%	21%	1%
Putrescine	66%	66%	<1%	22%	7%
Epinephrine	<1%	<1%	<1%	7%	92%
Lysine	2%	2%	<1%	6%	85%

*This data was cited from Ref.1

Application data


Polyamine imaging in both cancer and non-cancer cells by PolyamineRED

Three cancer cell lines (MCF7, MDA-MB-231 and SK-BR-3) and two non-cancer cells (MCF10A and human lymphocyte) were treated with 30 μM of PolyamineRED for 10 min. After incubation, cells were washed three times by PBS, followed by DAPI staining and formalin fixation. Images were obtained at Ex/Em=560 nm/585 nm for TAMRA and at Ex/Em=358 nm/461 nm for DAPI. TAMRA fluorescence was detected in cancer cell lines. On the other hand, incubation with non-cancer cell lines showed little fluorescence.



Evaluation of intracellular distribution of polyamines in MDA-MB-231 cancer cell lines

MDA-MB-231 cells were treated with 30 μ M of PolyamineRED for 10 min. After incubation, cells were washed three times by PBS, followed by DAPI-staining and formalin fixation. Images were obtained at Ex/Em=560 nm/585 nm for TAMRA and at Ex/Em=358 nm/461 nm for DAPI. Major TAMRA fluorescent signal was detected from nuclear. This indicates polyamines in MDA-MB-231 are mainly localized in nuclear.



Reference

1. K. Vong, K. Tsubokura, Y. Nakao, T. Tanei, S. Noguchi, S. Kitazume, N. Taniguchi, K. Tanaka, *Chem. Commun.*, **52**, 8403 (2017). Cancer cell targeting driven by selective polyamine reactivity with glycine propargyl esters.

Related product

AcroleinRED <Cell-based Acrolein Detection Reagent>


Acrolein is one of the most toxic oxidative stress marker and AcroleinRED is the world first cell-based acrolein detection reagent. As polyamines are one of the major source of acrolein, AcroleinRED and PolyamineRED are good set for oxidative stress research.

Catalog No. FDV-0022

Size 0.5 mg

Features

- Recommended Ex/Em: 560 nm/585 nm
- Easy and quick protocol
- Enable to semi-quantify intracellular acrolein



LipiRADICAL Green <Lipid Radical Detection Reagent>


LipiRADICAL Green is a specific fluorescent dye for lipid-derived radicals which are the most upstream factor of lipid peroxidation (LPO). LipiRADICAL Green can be applied into both *in vitro* assay and cell-based assay to monitor lipid radical productions.

Catalog No. FDV-0042

Size 0.1 mg

Features

- Recommended Ex/Em:~480 nm / 520 nm
- Enable to detect very unstable lipid-derived radicals
- Compatible with *in vitro* assay and in cell-based assay
- An innovative reagent for comprehensive identification of lipid-derived radicals by lipidomics



DNs-Rh <Cell-based GST Activity Assay Reagent >


Glutathione *S*-Transferases (GSTs) are major detoxication enzyme family which neutrize LPO-derived toxic aldehydes. DN_s-Rh is a novel fluorescent probe for monitoring wide GST members' activity both in cell and *in vitro*. DN_s-Rh releases green fluorophore rhodamine 110 upon GST activities. This probe has cell-permeability and can detect intracellular GST activity.

Catalog No. FDV-0031

Size 0.1 μmol

Features

- Easy and quick protocol
- Broad specificity for various GST family members
- Ex/Em: 496 nm/520 nm (commercial FITC filters are available)





E-mail Newsletter Sign Up

Our e-mail newsletter provides cutting edge information related to life science

Japanese



https://www.funakoshi.co.jp/users/register_user?fnads

English



<https://www.funakoshi.co.jp/exports/registration?ds>

 **funakoshi**
FRONTIERS IN LIFE SCIENCE

URL: <http://funakoshi.co.jp>
9-7 Hongo 2-Chome, Bunkyo-ku, Tokyo 113-0033