

Abstract

Distinct axonal and somatodendritic microenvironments can be established within compartmentalized XonaChips® by fluidically isolating one compartment from the other. This TechNote details this procedure using the dye, Alexa Fluor 488 hydrazide.

Introduction

Xona's compartmentalized devices and chips provide the unique ability to create distinct fluidic microenvironments¹. This TechNote details the procedure for establishing fluidic isolation between the main compartments of XonaChips®. For more information on cyclic olefin copolymer XonaChips®^{2,3}, see the [Introducing XonaChips®](#) technote.

Procedure

1. Prepare XonaChips® and culture neurons according to the appropriate protocol (refer to <http://www.xonamicrofluidics.com/protocols/>).
2. Remove 20 µL from the lower left well of the axonal compartment and place into the upper right well of the somatic compartment. Wait 2 min for flow within each channel to equilibrate.
3. Remove 50 µL of media from the axonal compartment. Add 0.3 µL of 1 mM Alexa Fluor 488 hydrazide to this media, mix via pipet and return back to the axonal compartment. The chip is ready for imaging.

NOTE: Other compounds of interest can be added. Adding a fluorescent dye with a similar molecular weight as the compound of interest is recommended in order to monitor


fluidic isolation over time.

Example Result

To demonstrate fluidic isolation, neurons were first dissociated from E18 rat hippocampus and cultured within two compartment XonaChips® (XC450) as described in the [XonaChip® protocol](#). Neurons were retrogradely labeled at 18 days *in-vitro*, by applying a G-deleted mCherry rabies virus²⁻⁶ to the axonal compartment. At 3 days post-infection, a low molecular weight dye- Alexa Fluor 488 hydrazide (pseudocolored magenta) was applied to the axonal compartment. No trace of fluorescence was detected within somatic compartment observed by confocal microscopy (Fig 1). Montage images of somatic and axonal compartments were acquired within one hour after adding the dye by laser scanning confocal microscope using a 30×/1.05 N.A. silicone oil (ne = 1.406) objective lens.

Conclusion

In summary, both XonaChips® and Xona's silicone devices provide fluidic isolation of axonal compartment from somatodendritic compartment. These fluidically isolated Xona devices provide compartmentalized application of drugs presenting high reproducibility and reliability in drug



screening studies.

If you are interested in testing XonaChips® contact us at info@xona.us

About Xona Microfluidics, Inc

Xona Microfluidics, Inc is a life sciences company based in Research Triangle Park, North Carolina. More information can be found at xonamicrofluidics.com.

References

1. Taylor, A. M. *et al.* A microfluidic culture platform for CNS axonal injury, regeneration and transport. *Nat Methods*. 2005; **2** (8), 599-605, doi:10.1038/nmeth777.
2. Nagendran, T., Poole, V., Harris, J. & Taylor, A. M. Use of Pre-Assembled Plastic Microfluidic Chips for Compartmentalizing Primary Murine Neurons. *J. Vis. Exp.* 2018; (141), e58421, doi:10.3791/58421.
3. Paranjape, S.R., Nagendran, T., Poole, V., Harris, J., and Taylor, A.M. Compartmentalization of Human Stem Cell-Derived Neurons within Pre-Assembled Plastic Microfluidic Chips. *J. Vis. Exp.* 2019; (147), e59250, doi:10.3791/59250.
4. Nagendran, T., Larsen, R.S., Bigler, R.L., Frost, S.B., Philpot, B.D., Nudo, R.J., and Taylor, A.M. Distal axotomy enhances presynaptic excitability onto injured pyramidal neurons via trans-synaptic signaling. *Nat Commun*. 2017 Sep 20; **8**(1):625. doi:10.1038/s41467-017-00652-y
5. Nagendran, T and Taylor, A.M. Unique axon-to-soma signaling pathways mediate dendritic spine loss and hyper-excitability post-axotomy. *Front. Cell. Neurosci.* 2019; Sep 24; **13**:431. doi: 10.3389/fncel.2019.00431.
6. Wickersham, I.R., *et al.* Monosynaptic restriction of transsynaptic tracing from single, genetically targeted neurons. *Neuron*. 2007 Mar 1; **53**(5):639-47.