


DynaMarker, MultiColor Stable Protein Ladders

Our "Protein MultiColor, Stable" can be stored at refrigerator (4°C) for 1 year !


Features

- **Stable at 4°C for 1 year.**
- The high stability keeps the sharpness of bands.
- M.W. range from approx. 17 kDa to 230 kDa.
- Ready-to-load (5 µL) prestain protein M.W. marker.
- Suitable for monitoring the mobility of each protein on SDS-PAGE gel and to assess blotting efficiency in western blotting without staining.



5 - 20% Polyacrylamide Gradient Gel

Comparison



SDS-PAGE after storing at appropriate temperature instructed by suppliers

- ★ BioDynamics Laboratory (BDL) "Protein Multicolor Stable" : Apply Volume : 5 µL
- 2. Company A : Apply Volume : 10 µL
- 3. Company B : Apply Volume : 10 µL
- 4. Company C : Apply Volume : 10 µL
- 5. Company D : Apply Volume : 5 µL

SDS-PAGE after heating at 37°C for 3 weeks

- ★ BioDynamics Laboratory (BDL) "Protein Multicolor Stable" : Apply Volume : 5 µL
- 2'. Company A : Apply Volume : 10 µL
- 3'. Company B : Apply Volume : 10 µL
- 4'. Company C : Apply Volume : 10 µL
- 5'. Company D : Apply Volume : 5 µL

➔ BioDynamics Laboratory's "Protein Multicolor Stable" can get very clear ladders even if heated at 37 °C for 3 weeks! Very Stable!

Supplier	BioDynamics Laboratory (BDL)	Company A	Company B	Company C	Company D
Apply Volume	5 µL	10 µL	10 µL	10 µL	5 µL
Storage Temp.	4°C	4°C	-20°C	4°C	-20°C
Storage Life	1 year	1 year	1 year	4 months	3 months


Product Information

[Manufacturer : BDL]

Product Name	Size	Code	Storage
Protein MultiColor, Stable, DynaMarker, Small	300 µL (60 loads)	DM660S	4°C
Protein MultiColor, Stable, DynaMarker	600 µL (120 loads)	DM660	4°C
Protein MultiColor, Stable, DynaMarker, Large	5 x 600 µL (600 loads)	DM660L	4°C

DynaMarker, Molecular Weight Protein Ladders

DynaMarker, Protein Eco (#DM610)



DynaMarker, Protein MultiColor III (#DM637)



Product Information

[Manufacturer : BDL]

Product Name	Size	Code	Storage
DynaMarker, Protein Eco	1 kit (120 loads)	DM610	-20°C
DynaMarker, Protein MultiColor III	2 x 300 µL (120 loads)	DM637	-20°C
DynaMarker, Protein MultiColor III, Large	10 x 300 µL (600 loads)	DM637L	-20°C

NOTE

※ Numbers after “#” represents product code.
 ※ All products here are research use only, not for diagnostic use.
 ※ Specs might be changed for improvement without notice.

※ Company name and product name are trademark or registered mark.
 ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor



funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,
 Tokyo 113-0033 JAPAN
 Phone : +81-3-5684-6296
 Fax : +81-3-5684-6297
 Email : export@funakoshi.co.jp