

Protein MultiColor Stable, Low Range


NEW

For more information : http://www.funakoshi.co.jp/exports_contents/80773

Pre-stained protein ladder for “Low” molecular weight proteins (1.7 to 47 kDa), stable in refrigerator (at 4°C) for 1 year !

Features

- Multicolored pre-stained natural proteins
- Stable for up to 1 year (at 4°C)
- Ready-to-load (5 µL / lane)
- Sharp bands
- M.W. range from approx. 1.7 to 47 kDa



- Ovalbumin (47 kDa)
- Carbonic anhydrase (31 kDa)
- Trypsin inhibitor (22 kDa)
- Lysozyme (17 kDa)
- Aprotinin (8.8 kDa)
- Insulin beta (3.8 kDa)
- Bacitracin (1.7 kDa)

Gel : 16 % Polyacrylamide (C3 %) /Tris-Tricine-SDS

Comparison

Supplier	BioDynamics Laboratory (BDL)	Company A	Company B	Company C	Company D
Apply Volume	5 µL	5 µL	5 µL	10 µL	10 µL
Storage Temp.	4°C	-20 °C	-20°C	-20°C	4°C / -20°C
Shelf Life	1 year	1 year	3 months	N/A	½ year / 1 year

Product Information

[Manufacturer : BDL]

Product Name	Code	Size	Storage
Protein MultiColor Stable, DynaMarker, Low Range	DM670	200 µL (40 loads)	4°C
	DM670L	3 x 200 µL (120 loads)	

Protein Multicolor Stable

For more information : http://www.funakoshi.co.jp/exports_contents/46111

Pre-stained protein ladder for “High” molecular weight proteins (17 to 230 kDa), stable in refrigerator (at 4°C) for 1 year !

Features

- Multicolored pre-stained natural proteins
- Stable for up to 1 year (at 4°C)
- Ready-to-load (5 µL / lane)
- Sharp bands
- M.W. range from approx. 17 to 230 kDa

Applications

- Real-time monitoring of protein mobility on SDS-PAGE gel
- Verification of western transfer efficiency on membranes without staining



- Myosin (230 kDa)
- β-Galactosidase (140 kDa)
- Phosphorylase-b (95 kDa)
- BSA (70 kDa)
- Ovalbumin (44 kDa)
- Carbonic anhydrase (32 kDa)
- Trypsin inhibitor (27 kDa)
- Lysozyme (17 kDa)

Gel : 5 ~ 20 % Polyacrylamide Gradient Gel

Comparison

Supplier	BioDynamics Laboratory (BDL)	Company A	Company B	Company C	Company D
Apply Volume	5 µL	10 µL	10 µL	10 µL	5 µL
Storage Temp.	4°C	4°C	-20°C	4°C	-20°C
Shelf Life	1 year	1 year	1 year	4 months	3 months

Product Information

[Manufacturer : BDL]

Product Name	Code	Size	Storage
Protein MultiColor Stable, DynaMarker	DM660S	300 µL (60 loads)	4°C
	DM660	600 µL (120 loads)	
	DM660L	5 x 600 µL (600 loads)	

NOTE

※ All products here are research use only, not for diagnostic use.
※ Specs might be changed for improvement without notice.

※ Company name and product name are trademark or registered mark.
※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor

Funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,
Tokyo 113-0033 JAPAN
Phone : +81-3-5684-6296
Fax : +81-3-5684-6297
Email : export@funakoshi.co.jp

BDL-1703-F08