



< Conventional procedure for DNA extraction >



- Complex and multiple procedure
- More than one hour

< Kaneka Easy DNA Extraction Kit version 2\* >



- Simple protocol
- 10 minute extraction


\*This product is for research use only.

< DNA extraction from knockout(KO) mouse and following PCR >

Extract DNA from mouse tails by using Kaneka Easy DNA Extraction Kit version 2\* following PCR.



Electrophoresis for Confirming DNA amplification



Kaneka DNA Easy Extraction Kit version 2 successfully extracted DNA from KO mouse tails.

