



< Kaneka Easy DNA Extraction Kit version 2* >



- Simple protocol
- 10 minute extraction

< Example of DNA extraction from shrimps >



< DNA extraction from shrimps and following PCR >


(Data obtained by our company)



1. Extract DNA from 20mg powdered shrimps by using Kaneka Easy DNA Extraction Kit version 2*.



Electrophoresis for Confirming DNA amplification



Kaneka DNA Easy Extraction Kit version 2 successfully extracted DNA from fish shrimps and Kaneka High-speed DNA Polymerase can amplify DNA from the extract.

* "Kaneka Easy DNA Extraction Kit version 2 is for research use only."