


# NIPRO Cell Gatch

**NIPRO Cell Gatch** is a protective exterior pouch for cryotubes, made of special plastic film, which remains durable and pliable in liquid nitrogen. NIPRO Cell Gatch protects from liquid nitrogen damage to cryotubes, eliminates cross-contamination within liquid nitrogen tank storage and during thawing. All cryotubes should be stored in this transparent protective pouch.



More Information :

[https://www.funakoshi.co.jp/exports\\_contents/80659](https://www.funakoshi.co.jp/exports_contents/80659)



**funakoshi**  
FRONTIERS IN LIFE SCIENCE

[www.funakoshi.co.jp/exports](http://www.funakoshi.co.jp/exports)  
Research Use Only



Our products make your cells happy!


# Cell Culture Essential Products

Funakoshi Co., Ltd. (Tokyo, Japan) offers many cell culture related products. These products will accelerate your research and helping to reveal new findings!

See more information at <https://www.funakoshi.co.jp/exports/> or contact us at [export@funakoshi.co.jp](mailto:export@funakoshi.co.jp)

# Disposable Hemocytometer (cell counting slide)

**Disposable Hemocytometer** is a cell counting slide with four chambers and same grid pattern as standard hemocytometer. This is made of high quality plastic and designed for rapid cell loading. No washing which saves time and avoid contaminations or unexpected expensive glass hemocytometer accidents .



More Information :

[https://www.funakoshi.co.jp/exports\\_contents/80019](https://www.funakoshi.co.jp/exports_contents/80019)



# NIPRO Culture Bag

**NIPRO Culture Bag** has excellent gas permeability with large surface area and is suitable for the culture of suspension cells such as lymphocytes and hybridomas in a closed sealed system with safer conditions without the need for rocking bag. Expansion culture is also easy by connecting two bags with special connecting tube.



1,000 mL and 350 mL size



Needleless support



Connecting Tube

NOTE

- ※ All products here are research use only, not for diagnostic use.
- ※ Specs might be changed for improvement without notice.
- ※ Numbers after “#” represents product code.
- ※ Company name and product name are trademark or registered mark.
- ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor

Funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,  
Tokyo 113-0033 JAPAN  
Phone : +81-3-5684-6296  
Fax : +81-3-5684-6297  
Email : [export@funakoshi.co.jp](mailto:export@funakoshi.co.jp)

FUN-2004-F19 (2020.04)

More Information :

[https://www.funakoshi.co.jp/exports\\_contents/80156](https://www.funakoshi.co.jp/exports_contents/80156)




# CryoScarless DMSO-free / Stemcell Keep

**CryoScarless DMSO-free** is a cell freezing medium, serum-free, 100% synthetic formulation for a wide variety of cells and cell lines.

**StemCell Keep** is the cell freezing medium for human ES/iPS cells.

Both these products are free from DMSO, Serum, Proteins and Xeno components for high recovery rate and low cytotoxicity cell storage conditions.

DMSO-free   Serum-free   Protein-free   Xeno-free



More Information :

[https://www.funakoshi.co.jp/exports\\_contents/46011](https://www.funakoshi.co.jp/exports_contents/46011)



# Accutase® / Accumax

**Accutase®** and **Accumax** is a natural enzyme mixture, non-mammalian for detachment of adherent cells, tissue or cell clumps dissociation. Both works gently, promptly and replaces trypsin and collagenase without neutralizing.

	Accutase®	Accumax
Enzyme Activity	Normal	Strong
Reaction Termination / Neutralizing	Not Required	Not Required
Animal and Bacterial Component	Not Contained	Not Contained
Phenol Red	Contained	Not Contained

More Information :

**Accutase®**

[https://www.funakoshi.co.jp/exports\\_contents/80531](https://www.funakoshi.co.jp/exports_contents/80531)

**Accumax**


[https://www.funakoshi.co.jp/exports\\_contents/80536](https://www.funakoshi.co.jp/exports_contents/80536)



# ThelioKeep

**ThelioKeep** is a preservation medium for dermal and nerve tissues in 4°C. This medium is free from DMSO, Serum, Protein and Xeno-components. ThelioKeep maintains structural strength and viability of these tissues.

DMSO-free   Serum-free   Protein-free   Xeno-free



**Tissue after preservation for 2 weeks at 4C.**

Human corneal endothelium tissue is preserved for 2 weeks at 4°C while maintaining structure and excellent cell viability.

More Information :



[https://www.funakoshi.co.jp/exports\\_contents/46140](https://www.funakoshi.co.jp/exports_contents/46140)



# iPGell / Smear Gell

**iPGell** is a cellular section preparation jellified agent for floating cells, 3D or cultured cells. It becomes an easy task to make a paraffin/frozen block and sectioning by using jellified samples.

**Smear Gell** is a jelly-like reagent for quickly mounting low number cell samples onto a specimen glass slide.



More Information :

**iPGell**

[https://www.funakoshi.co.jp/exports\\_contents/62940](https://www.funakoshi.co.jp/exports_contents/62940)

**Smear Gell**

[https://www.funakoshi.co.jp/exports\\_contents/80449](https://www.funakoshi.co.jp/exports_contents/80449)

