

PRODUCT INFORMATION

For Research Use Only

Product Name : Sonicated Salmon Sperm DNA

Code No. : F012, F013

Size : 1 ml / tube F012 1 ml
 F013 5 ml (1ml × 5)

Concentration: 10 mg/ml

Description: Salmon Sperm DNA was shared by sonication.
DNA is dissolved in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

Quality Control: Intrinsic DNase activity was not detected in the product,
by the use of pUC18 plasmid DNA as a control.

Preparation: Heat at 95°C for 3 minutes and cool on ice before use.

Storage Condition: - 20 °C

Related Products:

DM253	DynaMarker® Prestain Marker for Small RNA Plus Prestained Molecular Weight Marker for RNA (20-100 bases)
DM260	DynaMarker® Prestain Marker for RNA High Prestained Molecular Weight Marker for RNA (200-8,000 bases)

Ver. 4.1