

Titansphere™ Phos-TiO for Large Volume Samples



Titansphere™ Phos-TiO for Large Volume Samples


Appropriate for Larger Scale/Volume Purifications

The same specialized bulk media used in our Phos-TiO Kit is available in 50 mg/3 mL and 100 mg/3 mL cartridges as an extension of the Phos-TiO product line.

Sample Loading Capacity

Description	Content	
Sample	Tyr(PO ₃ H ₂) - Angiotensin II	
Spin Tip Sorbent Mass/Tip Volume	50 mg/3 mL	100 mg/ 3 mL
Sample Loading Capacity	50 µg	100 µg

Typical Enrichment Protocol using Titansphere Phos-TiO for Large Volume Samples



Phos-TiO for Large Volume Samples Cartridges are intended for use with a desktop or other centrifuge. While some of the versions of Phos-TiO resemble pipette tips or SPE cartridges, these products are not intended for use with pipettes or SPE vacuum manifolds; the cartridge internal configuration and particle size of the TiO bulk media requires centrifugal elution of all solutions.

Ordering Information

Titansphere™ Phos-TiO for Large Volume Samples

Description	Sorbent Mass/Tip Volume	Qty	Cat. No.
Titansphere Phos-TiO	50 mg/3 mL	25 pcs	5010-21290
for Large Volume Samples	100 mg/3 mL	25 pcs	5010-21291