

Lumiprobe Limited

8:00 - 17:00 HKT

Suite 12, 3/F, Great Eagle Centre, 23 Harbour Road, Wan Chai

Hong Kong

Phone: +852-8009-38633

Mob.: +852-5929-0488 (for calls from Hong Kong),

+86-147-14316277 (for calls from China)

Email: hk@lumiprobe.com

Cyanine 7.5 NHS ester


http://www.lumiprobe.com/p/cy75-nhs-ester

Cyanine 7.5 is a NIR dye with very long-wave emission. This product is amine-reactive derivative of this fluorophore.


Near infrared radiation readily penetrates tissues and can be used for imaging applications in vivo.

This product can be used for the labeling of various biomolecules containing amine groups, such as proteins and peptides, to track their distribution in organism in vivo by NIR imaging.

Increased rigidity of central polymethyne moiety allows to increase quantum yield by 20% compared to parent structure of Cy7.5® NHS.



Cy7.5 NHS ester structure



Cy7.5 absorbance and emission spectra

General properties

Appearance: green powder

Mass spec M+ increment: 630.4
Molecular weight: 833.76

CAS number: 2708152-94-5 Molecular formula: $C_{49}H_{52}N_3BF_4O_4$

IUPAC name: 1H-Benz[e]indolium, 2-[2-[3-[2-(1,3-dihydro-1,1,3-trimethyl-2H-benz[e]indol-2-

ylidene)ethylidene]-1-cyclohexen-1-yl]ethenyl]-3-[6-[(2,5-dioxo-1-pyrrolidinyl)oxy]-6-

oxohexyl]-1,1-dimethyl-

Solubility: soluble in organic solvents (DMSO, DMF, dichloromethane), low solubility in water

Quality control: NMR ¹H, HPLC-MS (95%)

Storage conditions: Storage: 12 months after receival at -20°C in the dark. Transportation: at room

temperature for up to 3 weeks. Avoid prolonged exposure to light. Desiccate.

Spectral properties

Excitation/absorption maximum, nm: 788 ϵ , L·mol⁻¹·cm⁻¹: 223000 Emission maximum, nm: 808 Fluorescence quantum yield: 0.10

Cy® is a trademark of GE Healthcare.