

Convenient and Safe! Disposable Homogenizer BioMasher SP

For more information : http://www.funakoshi.co.jp/exports_contents/80603

BioMasher SP is a convenient disposable homogenizer. This is suitable for extraction of DNA, RNA and protein from animal, insect and plant tissues.




Features

- High homogenization efficiency
 - By texture on both pestle and inside of tube.
- High collection amount
 - Enlarged type of BioMasher II.
- Heat can be controlled
- No centrifuge
- PP : Autoclavable / POM : non-Autoclavable
- EOG-Sterilized or non-sterilized product are available.
- Much higher homogenization efficiency by using PowerMasher II (Pestle Motor).

Handling procedure

1. Put samples and buffer into tube and insert pestle.
2. Grind samples by pestle by hand or Power Masher II.



PowerMasher II

PowerMasher II is a pestle motor, optimized for BioMasher II, III, V and SP series.

Features

- Increase homogenization efficiency dramatically.
- Works on two AA batteries (Accessories).



Check a movie at
http://www.funakoshi.co.jp/exports_contents/80603

Product Information

■ BioMasher SP

[Manufacturer : NIP]

Product Name	Size	Code
BioMasher SP (non-sterilized)	40 sets	893161
	100 sets	893162
BioMasher SP (EOG-sterilized)	40 sets	893163
	100 sets	893164

■ PowerMasher II (Pestle motor)

[Manufacturer : NIP]

Product Name	Size	Code
PowerMasher II	1 piece	891300

NOTE

- * All products on here are for research use only, not for diagnostic use.
- * Products might be changed for improvement without notice.

- * All brand names and product names are registered trademarks of their respective companies.

- * Please contact your local distributors for orders, quote request and any inquiry.

Your Local Distributor

funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,
 Tokyo 113-0033, JAPAN
 Phone: +81-3-5684-6296
 Fax: +81-3-5684-6297
 Email: export@funakoshi.co.jp